



# Instruction Manual ibidi Heating System, Universal Fit

Version 2.2



10918 ibidi Heating System, Universal Fit for 1 Chamber







# Contact

# ibidi GmbH

Lochhamer Schlag 11 82166 Gräfelfing Germany

Phone: +49 89 / 520 46 17 - 0 Fax: +49 89 / 520 46 17 - 59

> E-mail: info@ibidi.de Internet: ibidi.com



# Contents

1	Prea	amble	6
	1.1	Introduction	6
	1.2	Safety Symbols	6
	1.3	Nomenclature	7
	1.4	Specifications	7
	1.5	Disclaimer	10
	1.6	Safety Considerations	10
	1.7	Regulatory Statement	11
	1.8	Limited Warranty	12
	1.9	Transporting the ibidi Heating System, Universal Fit	12
	1.10	Repairing the ibidi Heating System, Universal Fit	13
	1.11	Waste Disposal – WEEE/RoHS Compliance Statement	13
		1.11.1 EU Directive WEEE	13
		1.11.2 EU Directive RoHS	13
2	Inte	nded Use	14
			11
3	Prin	nciple	14
3 4		ipment	
			14
	Equ	ipment	14 15
	<b>Equ</b> 4.1	<b>ipment</b> Components of the ibidi Heating System, Universal Fit	<b>14</b> <b>15</b> 15
	<b>Equ</b> 4.1 4.2	ipment         Components of the ibidi Heating System, Universal Fit         Combination Options	<ul><li>14</li><li>15</li><li>16</li></ul>
	<b>Equ</b> 4.1 4.2 4.3	ipment         Components of the ibidi Heating System, Universal Fit         Combination Options         Temperature Controller	<ol> <li>14</li> <li>15</li> <li>16</li> <li>17</li> </ol>
	Equ 4.1 4.2 4.3 4.4	ipment         Components of the ibidi Heating System, Universal Fit         Combination Options         Temperature Controller         Heated Lid	<ol> <li>14</li> <li>15</li> <li>15</li> <li>16</li> <li>17</li> <li>18</li> </ol>
	Equ 4.1 4.2 4.3 4.4 4.5	ipment         Components of the ibidi Heating System, Universal Fit         Combination Options         Temperature Controller         Heated Lid         Heated Plate	<ol> <li>14</li> <li>15</li> <li>16</li> <li>17</li> <li>18</li> <li>19</li> </ol>



5	Ope	eration	22
	5.1	Installation and Connection of the Parts	22
	5.2	Start the ibidi Heating System, Universal Fit	23
	5.3	Setting Parameters in the Front Display	23
		5.3.1 Run/Stop Mode	24
		5.3.2 Incubation Parameters	24
		5.3.3 Setup Menu	24
	5.4	Temperature Adjustment in the Sample	26
	5.5	Sample Preparation	27
6	Incu	aControl Software	28
6 7		aControl Software	28 29
-			
-	Mai	intenance	29
-	<b>Mai</b> 7.1 7.2	intenance Disinfection and Cleaning	<b>29</b> 29
7	<b>Mai</b> 7.1 7.2	intenance Disinfection and Cleaning Influence of Ambient Temperature and Ventilation	<b>29</b> 29 29
7	<b>Mai</b> 7.1 7.2 <b>Trou</b>	intenance Disinfection and Cleaning Influence of Ambient Temperature and Ventilation	<ul><li>29</li><li>29</li><li>29</li><li>30</li></ul>



# 1 Preamble

#### 1.1 Introduction

This manual is your guide to using the ibidi Heating System, Universal Fit for cell culture experiments on an optical microscope. It instructs first-time users how to use the instrument, and serves as a reference for experienced users.

Before using the ibidi Heating System, Universal Fit, please read this instruction manual carefully and make sure that the contents are fully understood. This manual should be easily accessible to the operator at all times during instrument operation. If this manual gets lost, order a replacement from www.ibidi.com.

To ensure safe operation, the ibidi Heating System, Universal Fit must only be operated with the supplied components and according to the instruction manual.

## 1.2 Safety Symbols

Note that the signal words **WARNING**, **CAUTION** and **NOTE** have specific meanings in this manual. Do not proceed beyond a signal word until you have performed the indicated actions.

- **WARNING!** A potentially hazardous situation which, if not avoided, could result in serious injury or even death. Warning messages in the text are displayed in a gray shaded box.
- **CAUTION** A potentially hazardous situation which, if not avoided, could result in minor or moderate injury. It is also used to alert against damaging the equipment or the instrument.
- **NOTE** Additional information to help achieve optimal instrument and assay performance.

Symbols on the product identification label and back panel of the device:



CE Marking: This symbol indicates the product's compliance with EU legislation.



This label is positioned on the back of the device and prompts you to read the manual before using the device.



Product disposal: The symbol indicates that this product must be recycled/disposed of separately from other household waste. See page 13 for details.



#### 1.3 Nomenclature



# 1.4 Specifications

Only operate the Heating System in the range of the specifications given below:

Table 1 – Specifications	of the ibidi Heating System, Universal Fit
--------------------------	--

Electrical Specifications Power Supply		
Protection class	Ι	
International protection marking	IP 20	
(IEC 60529)		
Overvoltage category	II	
External power supply	AC 100-240 V, 50/60 Hz, 2 A	
Input line voltage Temperature	DC 24 V, 6.67 A, 160 W	
Controller		
Output voltage to channel 1	DC 10 V, max. 5 A	
(Heated Lid), Universal Fit		
Output voltage to channel 2	DC 12 V, max. 5 A	
(Heated Plate), Universal Fit, for 1		
Chamber		

Operating and Storage Condi		
Operating site	Indoor use only	
Operating temperature	18-30°C/64-86°F (min 5°C/9°F less than set temperature)	
Operating humidity	max. 80% relative humidity (RH)	
Operating altitude	max. 2000 m (atmospheric pressure 800-1060 hPa/11.6-15.4 psi)	
Storage conditions	-5-50°C/23-122°F, humidity <60% relative humidity (RH)	

Outer Dimensions and Characteristics of the Components		
Temperature Controller90 mm × 170 mm × 230 mm		
	Weight: 1720 g/3.8 lbs	
Heated Lid	19 mm × 85.5 mm × 127.5 mm (134.5 mm with cover ridge) Length of cable: 1.5 m	



	Table 1 – (continued)
Heated Plate Heated Plate assembled with Heated Lid	Connector to Gas Incubation: Female Luer Lock 12 mm × 85.5 mm × 127.5 mm Observation area: 40 mm × 82 mm Length of cable: 1.5 m Height: 25.5 mm
Heating Inserts all Heating Insert 35 mm Dish high	Weight: 330 g/0.46 lbs 47 mm × 97 mm Height without blank holder: 8 mm Height with blank holder: 16-18 mm Observation area: ∅22 mm Weight insert: 50 g/0.11 lbs Weight blank holder: 72 g/0.16 lbs
Heating Insert 35 mm Dish low	Height without blank holder: 8 mm Height with blank holder: 13-15 mm Observation area: Ø22 mm Weight insert: 50 g/0.11 lbs Weight blank holder: 47 g/0.11 lbs
Heating Insert Slides	Height without blank holder: 8 mm Height with blank holder: 13-15 mm Observation area: 49 mm × 22 mm Weight insert: 42.5 g/0.94 lbs Weight blank holder: 38.5 g/0.85 lbs
USB cable Power supply cable	Length: 1.8 m Length: 2.0 m (power supply to wall) Length: 1.2 m (power supply to device)

Temperature Control Range	
Heated Lid	Ambient temperature (min. $18^{\circ}$ C) to $+45^{\circ}$ C Accuracy: $\pm 0.2^{\circ}$ C (at sensor location)
Heated Plate	$\pm 5^{\circ}$ C (entire heated glass) Ambient temperature (min. 18°C) to +45°C Accuracy: $\pm 0.2^{\circ}$ C (at sensor location) $\pm 0.5^{\circ}$ C (entire heated plate)

## **Recommended Temperatures before Adjustment**

See Section 5.4

#### **USB** Interface

Connector type	USB 2.0 Connector Type B
Recommended USB cable	Tripp Lite UR022-006 (shielded)
Driver	FTDI VCP driver



## Table 1 – (continued)

Microscope Requirements			
Working distance condenser	$\geq$ 26 mm		
Stage holder	Holder for standard multi-well plates		
<b>Optical Properties Glass Lid</b>	Optical Properties Glass Lid		
Glass	Selected HQ Floatglass		
Thickness of the glass plate	1.1 mm		
Refractive index glass	$n_D = 1.520 (588 \text{ nm})$		
ITO coating	Thickness: 100 nm		
Passivation layer	Thickness: 20-25 nm		
Refractive index ITO with passiva-	1.95		
tion layer			



#### 1.5 Disclaimer

- ibidi shall not be held liable, either directly or indirectly, for any damage incurred as a result of product use.
- The contents of this manual are subject to change without notice for product improvement.
- This manual is considered complete and accurate at publication.
- This manual does not guarantee the validity of any patent rights or other rights.
- If an ibidi software program doesn't function properly, this may be caused by a conflict from another program operating on the computer. In this case, take corrective action by uninstalling the conflicting product(s).
- ibidi is a registered trademark of ibidi GmbH in Germany and other countries.

## 1.6 Safety Considerations

#### WARNING!

- Only operate the ibidi Heating System, Universal Fit with the supplied components.
- Only use the cables and plugs delivered with the system. The power plug of the control unit must be inserted in an outlet with a ground (earth) contact.
- Do not replace detachable power cables by power cables with inadequate specifications. By violating these instructions you risk electric shock and fire.
- Only use extension cables that have a protective ground wire.
- Do not operate the ibidi Heating System, Universal Fit under conditions that pose a risk of explosion, implosion, or the release of gases. Only operate the ibidi Heating System, Universal Fit with aqueous solutions.
- Do not operate a damaged ibidi Heating System, Universal Fit. If the housing seems damaged or something is rattling inside the controller, contact the ibidi service hotline for repair.
- Some accessible parts of the Heated Plate and Heated Lid can reach temperatures up to 55°C. Avoid touching the temperature-controlled parts of the system when you have set the Temperature Controller to high temperatures.

#### CAUTION

- Ensure that the external power supply is easily accessible. The ibidi Heating System, Universal Fit must be installed in a manner such that none of its components hinders access to the external power supply.
- Immediately replace damaged cords, plugs, or cables to avoid risk of personal injury or damage to the instrument.



- Only ibidi technical staff and technical staff instructed by ibidi are permitted to open and service the ibidi Heating System, Universal Fit.
- The external power supply should not be brought into contact with moisture. If the housing is damaged, the external power supply should not be used.
- Avoid strong magnetic fields and sources of high frequency. The ibidi Heating System, Universal Fit might not function properly when located near a strong magnetic field or high frequency source.
- Avoid vibrations from vacuum pumps, centrifuges, electric motors, processing equipment, and machine tools.
- Avoid dust and corrosive gas. Do not install the ibidi Heating System, Universal Fit where it could be exposed to high levels of dust or to outside air or ventilation outlets.
- Install the ibidi Heating System, Universal Fit in a horizontal and stable position, such as a table, bench, or desk upon which the instrument is installed.
- Install the ibidi Heating System, Universal Fit in a location that enables easy access for maintenance.
- Do not place heavy objects on the instrument.
- The heated glass plates of the incubation chamber can break on mechanical impact. If so, the glass shards can lead to injuries if handled.
- Be aware that when switched on, a 10 V DC voltage is applied to the underside of the glass on the Heated Lid. Do not touch the underside or put it in contact with anything conductive. This could cause a short circuit that may damage the Temperature Controller and/or the Heating Devices.

# 1.7 Regulatory Statement

The ibidi Heating System, Universal Fit has been designed, produced and tested in compliance with the European standard DIN EN 61010-1 (IEC 61010-1, "Safety requirements for electrical equipment for measurement, control and laboratory use"). Furthermore it meets the IEC 61326-1 ("Electrical equipment for measurement, control and laboratory use - EMC requirements") and CISPR 11 ("International Standard for electromagnetic emissions (disturbances) from Industrial, Scientific and Medical (ISM) Equipment") standards .

The device carries the CE mark.

The ibidi Heating System, Universal Fit meets the Low Voltage Directive 2014/35/EU and the EMC Directive 2014/30/EC.



#### 1.8 Limited Warranty

Products manufactured by ibidi, unless otherwise specified, are warrantied for a period of one year from the date of shipment to be free of defects in materials and workmanship. If any defects in the product are found during this warranty period, ibidi will repair or replace the defective part(s) or product free of charge.

This warranty does not apply to defects resulting from the following:

- 1. Improper or inadequate installation.
- 2. Improper or inadequate operation, maintenance, adjustment, or calibration.
- 3. Unauthorized modification or misuse.
- 4. Use of unauthorized tubing or fluidic connectors.
- 5. Use of consumables, disposables, and parts not supplied by an authorized ibidi distributor.
- 6. Corrosion due to the use of improper solvents, samples, or due to surrounding gases.
- 7. Accidents beyond ibidi's control, including natural disasters.

This warranty does not cover consumables, such as cell culture chambers and dishes, tubes, fluidic connectors, reagents etc.

The warranty for all parts supplied and repairs provided under this warranty expires on the warranty expiration date of the original product.

## 1.9 Transporting the ibidi Heating System, Universal Fit

The weight of the Temperature Controller is approx. 1.7 kg/3.8 lbs. Moving the Temperature Controller during operation can pose a risk of personal injury or damage to the instrument.

For transport, switch off the Temperature Controller and then disconnect the heated components from the controller. Carry the devices carefully and avoid mechanical shocks.

#### WARNING!

Hot surface (max. 55°C)! Do not touch Heated Lid and Heated Plate when hot. Always disconnect the instrument from the power supply before transport and leave the instrument to cool down for approx. 5 minutes.



#### 1.10 Repairing the ibidi Heating System, Universal Fit

For inquiries concerning repair service, contact the ibidi service personnel and provide the model name and serial number of your system.

ibidi GmbH Service Hotline: service@ibidi.com

## CAUTION

Do not try to repair the ibidi Heating System, Universal Fit by yourself. Disassembly of the ibidi Heating System, Universal Fit is not allowed. Disassembly poses a risk of personal injury or damage to the devices. Contact ibidi service personnel if there is a need to disassemble a device.

#### 1.11 Waste Disposal – WEEE/RoHS Compliance Statement

The European Union (EU) has enacted two directives, the first on product recycling (Waste Electrical and Electronic Equipment, WEEE) and the second on limiting the use of certain substances (Restriction on the use of Hazardous Substances, RoHS).

#### 1.11.1 EU Directive WEEE

The ibidi Heating System, Universal Fit must be disposed of in compliance with the WEEE Directive 2012/19/EC.



This symbol on the product is in accordance with the European Union's Waste Electrical and Electronic Equipment (WEEE) Directive. The symbol indicates that this product must be recycled/disposed of separately from other household waste. It is the end user's responsibility to dispose of this product by taking it to a designated WEEE collection facility for the proper collection and recycling of the waste equipment. The separate collection and recycling of waste equipment will help to conserve natural resources and protect human health and the environment. For more information about recycling, please contact your local environmental office, an electrical/electronic waste disposal company or distributor where you purchased the product.

#### 1.11.2 EU Directive RoHS

The ibidi Heating System, Universal Fit meets the requirements set forth in the RoHS Directive 2011/65/EU and the additional terms as specified in RoHS Directive 2015/863/EU.



# 2 Intended Use

The ibidi Heating System, Universal Fit is a stage top incubator for live cell imaging that fits onto inverted microscope stages with a universal mounting frame for multiwell plates. The geometry of the chamber, with a heated glass plate above the sample, provides a platform for keeping microscopy slides and dishes at a constant temperature during an experiment on a microscope.


For full control of the incubation environment, the ibidi Heating System, Universal Fit can be combined with the ibidi Gas Incubation System to regulate the  $CO_2$  and  $O_2$  concentrations as well as the humidity.


The modular concept of the Heating Inserts enables the use of a variety of microscopy labware ranging from 35 mm dishes to regular slides.


# 3 Principle


## **Physiological Conditions in Live Cell Imaging**

The Heated Lid and Heated Plate are designed to keep cells on-stage at 37°C in microscopy slides and dishes (Figure 1). To achieve this, the glass top of the Heated Lid and the Heated Plate are actively heated by the Temperature Controller. The Heated Lid prevents condensation effects inside the incubation system (Figure 2).












# 4 Equipment

This section provides a brief overview over all parts of the ibidi Heating System, Universal Fit including a description of the characteristics.

## 4.1 Components of the ibidi Heating System, Universal Fit

The components of the ibidi Heating System, Universal Fit are listed below.





Description	Drawing
Heating Insert µ–Slide*	
Heating Insert LabTek*	
Temperature Adjustment Set: 1 temperature sensor 1 perforated ibidi μ–Slide 8 Well 1 perforated ibidi μ–Dish <sup>35mm, high</sup> 1 perforated ibidi μ–Dish <sup>35mm, low</sup>	
USB cable to connect the Temperature Controller with a computer	
Country specific power cord to connect the external power supply to the wall socket	
External power supply for the Temperature Controller	
USB flash drive with IncuControl software	L C Line Andrew C

\*One Heating Insert of your choice is delivered with the ibidi Heating System, Universal Fit. If needed, more Heating Inserts can be ordered separately.

# 4.2 Combination Options

The parts of the ibidi Heating System, Universal Fit are combined as shown in Figure 3. The Heated Lid and Heated Plate fit on an inverted microscope stage equipped with a universal mounting frame for multiwell plates. They are connected to the Temperature Controller.

The ibidi Heating System, Universal Fit can be combined with the Gas Mixer unit, that provides  $CO_2$  and  $O_2$  (optional) control and a defined humidified atmosphere. Detailed information on the Gas Mixer unit is given in the Gas Mixer instructions.

To hold several geometries of Slides and Dishes, the Heated Plate can be equipped with different Inserts (Section 4.6): Insert for  $\mu$ -Dish <sup>35mm, high</sup>, Insert for  $\mu$ -Dish <sup>35mm, low</sup>, Insert for  $\mu$ -Slides, and Insert for LabTek<sup>TM</sup> chambered coverglass.

For parallelization of experiments, the Heated Plate for 4  $\mu$ -Slides provides a platform to observe 4 Slides in parallel. Detailed information is given in the instructions of the ibidi Heating System, Universal Fit for 4  $\mu$ -Slides (#10927).

#### ibidi Heating System, Universal Fit






Figure 3 – Overview of the parts of the ibidi Heating System, Universal Fit with options to combine

## 4.3 Temperature Controller

The Temperature Controller is designed to control the different heated components (four channels available).

The front display shows the set values (S) and the current values (I). The settings can be adjusted via the control buttons (Section 5.3).

If you wish to control the Temperature Controller via PC, use the IncuControl software (Section 6).






Figure 4 – Front view of the Temperature Controller.

The LEDs on the front indicate the status of the channels, connection to the power supply, and USB connection (Figure 4).

Control LEDs Channel 1-4		
LED off Channel inactive		
LED on	Channel active	
LED fast blinking	Channel error	

All plugs for the electrical connections are integrated into the rear of the Temperature Controller (Figure 5). The setup of the connections is shown in Section 11.




Figure 5 – Rear view of the Temperature Controller.


# 4.4 Heated Lid

The Heated Lid provides excellent optical quality, and also allows for the use of all standard microscopy techniques, including differential interference contrast (DIC). Due to the height of the lid, we recommend using condensers with a working distance of  $\geq 26$  mm.



The upper glass part of the Heated Lid is heated. The electrical cable is connected to the Temperature Controller. The Heated Lid fits exactly and securely onto the Heated Plate.

For gas incubation, the Heated Lid is equipped with inlets for the gas flow and the humidity sensor. Detailed information on the Gas Incubation System is given in its instruction manual.



# 4.5 Heated Plate

The Heated Plate provides the base for the ibidi Heating System, Universal Fit. The lower part of the Heated Plate fits into any microscope stage with a universal mounting frame for multiwell plates. The electrical cable (1.5 m) is connected to the Temperature Controller.




Figure 7 – Heated Plate, top view



## 4.6 Heating Inserts

The Heating Insert holds the sample and fits into the Heated Plate. The Insert is held in place by strong neodymium magnets. Use the handles to place/remove the insert into/from the Plate of the ibidi Heating System, Universal Fit.

The Heating Inserts have a lower part (insert) that fit into the Heated Plate, and an upper part (holder) to press the sample down. The two-part inserts have two functions that both rely on the force of the integrated magnets. The first function is to hold the Slides or Dishes firmly in position, so as to avoid displacement during microscope stage movements. The second is to create a tight contact between the Slides or Dishes and the metallic insert, to maximize the heat transfer and ensure stable heating of the sample. The Heated Plate passively heats the inserts.

#### WARNING!

Heated Plate, Heated Lid, and Heating Inserts all contain strong neodymium magnets! Please contact us for a non-magnetic system, if permanent magnetic fields are detrimental to your experiment.

#### 4.6.1 Heating Insert µ–Slide

The Heating Insert  $\mu$ -Slide fits 75 mm × 25 mm  $\mu$ -Slides (e.g. all ibidi  $\mu$ -Slides).




Figure 8 – Heating Insert µ–Slide

## 4.6.2 Heating Insert $\mu$ -Dish <sup>35mm, high</sup> and Heating Insert $\mu$ -Dish <sup>35mm, low</sup>

The Heating Insert  $\mu$ –Dish <sup>35mm, high</sup> fits only ibidi  $\mu$ –Dish <sup>35mm, high</sup>, the Heating Insert  $\mu$ –Dish <sup>35mm, low</sup> fits only ibidi  $\mu$ –Dish <sup>35mm, low</sup>

Center the  $\mu\text{-Dish}$  in the insert before you put the holder on.






Figure 9 – Heating Insert  $\mu$ –Dish

To make sure that the  $\mu$ -Dishes are properly positioned, the holder of the insert must apply a diagonal force in order to snap the dishes into a fixed position. This will leave an inclined gap between the holder and the insert, as shown in Figure 10, but this is intentional and helps to keep the dishes in place.



Figure 10 – Schematic Drawing of  $\mu$ -Dish in the Heating Insert. To ensure a tight fit, the lid of the Insert is tilted when closed.



# 5 Operation

Before starting an experiment, check that the ibidi Heating System, Universal Fit fits on your microscope stage, and that your cell culture vessels are compatible with the Heating Insert. Connect all parts (Section 5.1) and perform a temperature adjustment as explained in Section 5.4.

# NOTE

The Temperature Controller only measures the temperature from the sensors in the Heated Lid and Heated Plate. The temperature in the sample must be adjusted for your specific setup. Follow the instructions in Section 5.4.

## 5.1 Installation and Connection of the Parts

The components of the ibidi Heating System, Universal Fit are connected as shown in Figure 11.




Figure 11 – Installation of the components of the ibidi Heating System, Universal Fit.

# NOTE

Before shipment, all controllers are run through an in-house calibration with the corresponding Heated Lid and Heated Plate. Only use the Temperature Controller with the corresponding Heated Lid and Heated Plate to ensure correct sensor calibration.

- 1. Place the Temperature Controller next to the microscope and connect the power supply.
- 2. Confirm that the power switch of the Temperature Controller is off.



- 3. Insert the Heated Plate into the opening of the microscope stage.
- 4. Insert an empty, unsterile Slide/Dish into the Heating Insert and put the assembly in the opening of the Heated Plate to close it during equilibration.
- 5. Place the Heated Lid onto the Heated Plate.
- 6. Plug the electrical cables of the Heated Lid (Channel 1) and Heated Plate (Channel 2) into the connectors on the back of the Temperature Controller.
- 7. Let the system equilibrate for at least 30 minutes.

To setup communication with the IncuControl Software, the USB cable must be connected between the Temperature Controller and the computer.

To adjust the sample temperature, use the temperature sensor in the Temperature Adjustment Set and plug it into the corresponding connector on the back of the Temperature Controller. The other end is placed in the sample (Section 5.4).

#### 5.2 Start the ibidi Heating System, Universal Fit

The Temperature Controller is switched on by the dip–switch on the back. Make sure the display shows "RUN" (see Section 4.3). The system immediately starts heating up the heated components.

#### WARNING!

Be aware that when the system is switched on, 10V DC voltage is applied to the underside of the glass plate. Do not touch the underside or contact it with anything conductive! This could cause a short circuit that may destroy the controller and/or the lid.

If you are not sure which temperature is set (e.g. when operating the system the first time), it is recommended to disconnect all heated devices (Heated Plate, Heated Lid, and other optional heated devices) and then switch the system on. By doing this, you ensure that you do not start the heating process with the wrong temperature settings. It is now possible to set the temperatures for the individual channels (Section 5.3). After this, you can re-connect the heated devices.

#### 5.3 Setting Parameters in the Front Display

All control parameters can be manually set on the controller using the buttons and the display on the front panel.




Figure 12 – Temperature Controller front display and set buttons.



The display shows the measured ('I' = instantaneous) and set ('S' = set) temperatures of all the channels.

The cursor position is indicated with square brackets ("[ ]"). You can move the cursor using the "left" and "right" buttons. If you want to select a parameter or a function, press the round button and the square brackets will change to angle brackets ("< >"). Now you are able to change the value with the "up" and "down" buttons. To confirm the changed value, you must press the round button once more.

Standardly, the channels are pre-defined with the following settings to match all possible combinations of heated components on the controller:

Channel No.	Short name	Components to connect	Product number
Channel 1	"Lid"	Heated Lid, Universal Fit	#10918, #10927
		or Heated Lid, Multi-Well Plates, K-Frame	#10929
Channel 2	"Plate"	Heated Plate, Universal Fit, for 1 Chamber	#10918
		or Heated Plate, Multi-Well Plates, K-Frame	#10929
Channel 3	"Glass"	Heated Glass Bottom, Multi-Well Plates, K-Frame	#10929
Channel 4	"Ch4"	Heated Plate, Universal Fit, for 4 µ-Slides	#10927, #10928

Table 4 - Channel assignment of the Temperature Controller

#### 5.3.1 Run/Stop Mode

Set the whole system to run or stop mode by manipulating the setting in the display's bottom line to "Run" (= system is running) or "Stop" (= system is not running).

	Lid	Plate	Glass
Ι	42.1	37.0	N.C.→
<b>S</b> [	42.0	] 37.0	42.0
	N.C.	Run	Setup

Figure 13 – Main display indicating "Run"

#### 5.3.2 Incubation Parameters

Set the parameters in the front display and wait for the temperature to equilibrate. It is recommended to start the system at least 30 minutes before inserting cells. Recommended values are given in table 5.

Before starting an experiment, an initial temperature adjustment must be performed (see Section 5.4).

#### 5.3.3 Setup Menu

The Setup menu offers the possibility to change the settings of the individual channels and the display. For standard operation no changes need to be made. Enter the setup menu by navigating to the "Setup" entry and pressing the round button. The mode dialog will open.



**Mode** Each channel of the Temperature Controller can be set to an "On" or "Off" state. For standard operation all channels must be turned on.

Setup		<mode></mode>
Channel	=	1
Mode	=	On
		Return

Figure 14 – Mode dialog

- 1. Press the round button once more to move the cursor to the channel number.
- 2. Select the respective channel and confirm with the round button.
- 3. Select "On" or "Off" and confirm with the round button.
- 4. To move on navigate to "Return" and press the round button.

**Alarms** In this dialog it is possible to set the alarm limits for divergent control parameters. Move from the mode menu to the alarm menu by pressing the "right" button.

The maximum and minimum limits of the alarm can be set for each individual channel. If the alarm is activated (which happens when the current value goes under/over the low/high limits), the current value (I) blinks, showing alternately the value and "low" or "high", respectively. To stop the display blinking, navigate to the blinking channel and press the round button once.

Setup		<alarms></alarms>
Channel	=	4
High	=	42°C
Low	=	36°C

Figure 15 – Alarm settings dialog

**Preferences** Set the brightness and contrast of the display in the preferences menu.

Setup <preferences></preferences>
Backlight = 65
Contrast = 15
Return

Figure 16 – Preferences dialog

Info Info about serial number and firmware version is shown in the info dialog.

Setup	<info></info>
2 SN: ibiTC2	0.01
FW 2.11.02	-001

Figure 17 – Info display



## 5.4 Temperature Adjustment in the Sample

The ibidi Heating System, Universal Fit controls the temperature of the actively heated components. The temperature in the sample results from the interaction of all heated parts. Therefore, the temperature at the position of the cells must be measured and adjustments should be made to the heated components to achieve the desired temperature. An adjustment is recommended before starting the first experiment. The temperature adjustment should be done for each chamber type, in order to control temperatures with an absolute accuracy of less than 1°C.

When starting the temperature adjustment in the sample in your specific setup, use the recommended temperature settings (for samples that should be 37°C):

Table 5 – Recommended temperature settings before optimization

Device	Temperature
Heated Lid (channel 1)	42°C
Heated Plate (channel 2)	38°C

To measure and calibrate the sample temperature, use the provided Temperature Adjustment Set. The temperature of the heated components has to be adjusted in small steps. Follow this procedure to adjust the sample temperature:

- 1. Set all your peripheral experimental parameters, such as room temperature, air conditioning, airflow, illumination, microscope settings, etc.
- 2. Connect the temperature sensor (thermocouple type K) to the plug on the Temperature Controller's back.
- 3. Fill the provided μ–Dish or μ–Slide (depending on what you want to use for your experiments) with water. Volumes should resemble the amount you need for your experiments.
- 4. Put the loose end of the thermocouple through the hole in the lid of the  $\mu$ -Dish/ $\mu$ -Slide and make sure it dips into the water.
- 5. Place the  $\mu$ -Dish/ $\mu$ -Slide in the Insert and put the holder on.
- 6. Place the Insert in the ibidi Heating System, Universal Fit and close the Heated Lid.
- 7. After 30 minutes, check the temperature of the external sensor on the display or in the Incu-Control software. If the sample temperature is still too low, raise the set values for the Heated Plate and Heated Lid for 0.5-1°C and wait until the sample temperature is stable again (at least 10-15 min).
- 8. When the sample temperature has reached the desired value, write down the set values for the Heated Plate and Heated Lid and use those settings for all upcoming experiments with the same probe holder.



## NOTE

The temperature of the Heated Lid must be set to at least some degrees warmer than the temperature of the Heated Plate!

The temperature adjustment must be repeated from time to time (at least once a year), especially if one of the following conditions has been changed:

- Room temperature
- Air conditioning
- Chamber type or objective lens
- Humidification and gas flow
- Use of an XL-Incubator
- Use of an Objective Heater

#### 5.5 Sample Preparation

Check that the cell cultureware you intend to use fits into the Insert of the ibidi Heating System, Universal Fit.

Prepare the cells according to your protocol and place the culture vessel in the Heating Insert. The Heating Insert can be removed from the ibidi Heating System, Universal Fit even with the holder in place. Then put the whole assembly in the opening of the Heated Plate and close the Heated Lid.

## NOTE

Let the temperature of the system equilibrate for a minimum of 30 min before you start your experiments.



# 6 IncuControl Software

The Temperature Controller has a USB interface for computer control and data logging. For this purpose, ibidi provides the IncuControl software that comes with the controller or can be downloaded from the ibidi website.

For more details, please refer to the IncuControl instructions.

DataLogViewer	About	Alarm settings	Offsets		
Power ON		Tem	oControl		<b>Whidi</b>
		Ch 2		Ch 4	cells in focus
Title	Lid 💌	Plate 💌	Lid 💌	Plate 💌	
Measured	39,9 °C	37,0 °C	0,0 °C	0,0 °C	
Set Value	<b>40,0 🚔</b> "C	<b>37,0 🚔</b> °C	<b>40,0 🚔</b> °C	37,0 🚔 °C	
	ON			ON	
		Record		Ther	mo S. 37 °C
Da 100-	ta Recording	1		ON Clear	]
					Ch1 📈
- <sup>08</sup>					Ch 2 📈 Ch 3 📈
arnte 60 -					Ch 4 📈
- 08 - 09 - 09 - 09					Ch 5 📈
<sup>₩</sup> 20-					
0-					
					05:00
	🔵 сом	0			

Figure 18 – Temperature control window in IncuControl.

Prohibitions on the use of ibidi software:

- Copying software for other than backup purposes
- Transfering or licensing of the right to use software to a third party
- Disclosure of confidential information regarding software
- Modification of software



# 7 Maintenance

#### 7.1 Disinfection and Cleaning

The heated components and the Temperature Controller do not have to be sterile, because the incubation chamber has no direct contact with the cells and the cell culture medium. If disinfection is necessary for some reason, we recommend using isopropanol (70%) or common lab disinfection solutions based on quaternary ammonium compounds (e.g. Barrycidal 36 or Pharmacidal).

All parts of the ibidi Heating System, Universal Fit can be cleaned from the outside. We recommend using ultrapure water for cleaning. Fingerprints on the Heated Lid can be removed using isopropanol (70%) or lens/eyeglass cleaning wipes.

## CAUTION

When cleaning the heated glass plate(s), be careful when wiping the inner surface of the lid not to damage the electro-conductive coating. Also take care not to damage the glass plate. The use of ethanol or other types of organic solvents may remove the instrument's paint.

To clean the heated components and/or the Temperature Controller switch off the Temperature Controller and disconnect all electrical cables. Leave the instrument to cool down for approx. 5 minutes.

#### WARNING!

Hot surface (max. 55°C)! Do not touch Heated Lid, Heated Plate, and Heated Glass Bottom when hot.

## 7.2 Influence of Ambient Temperature and Ventilation

The ambient temperature affects the conditions inside the incubation chamber. Devices, such as computers and camera controllers, can significantly heat up small rooms. In this case, we recommend equilibrating the room temperature to the typical experimental conditions at least 2-3 hours before starting the experiment.

Ventilation can enhance the effect of temperature and humidity changes in the vicinity of the incubation chamber. In a case where the airflow (e.g., air conditioning) cannot be stopped, we recommend protecting the microscope as much as possible.

The use of an XL-Incubator and/or an objective heater minimizes those effects and helps significantly to stabilize surrounding conditions.



# 8 Troubleshooting

## 8.1 Focus not Stable

Focus drift is detrimental for most microscopy experiments, especially long duration time-lapse experiments. Focus stability is mainly influenced by mechanical changes and temperature variations. Follow these recommendations to keep your cells in focus:

- Switch on all components (e.g., heating, gas incubation, computer, or other equipment) at least 60 minutes before starting the experiment.
- After you put the sample onto the microscope, wait 30 minutes before starting a time-lapse experiment to achieve temperature and immersion oil equilibration<sup>1</sup>.
- Keep the room temperature as stable as possible. Air conditioning should either be working continuously or switched off.
- Do not change the temperature during the experiments. Avoid door/window openings, as this could rapidly change the temperature.
- Eliminate all sources of mechanical vibrations. Use a damped table for your microscope.

## 8.2 Evaporation Is too High

Depending on the incubating conditions, small volumes might evaporate quickly, especially during long-term experiments. If you have an actively controlled humidifying device (e.g. ibidi Gas Incubation System), increase the set value for relative humidity. Additionally, we suggest using silicone oil (e.g. Anti-Evaporation Oil, ibidi, 50051) to decrease evaporation.

Covering the medium with sterile silicone oil prevents all evaporation effects and is compatible with cell culture. Please do not use mineral oil, as this can be harmful to your cultureware.

Equilibrate oil and medium inside the incubator overnight. This step helps to avoid the formation of air bubbles, and pre-warms the solutions to 37°C. Afterwards, fill your slide with cells and medium. Cover the medium's surface with an appropriate amount of silicone oil. Don't drip the oil directly onto the surface, but let it run down the edges of the cell culture vessel. Details about avoiding evaporation are given on the ibidi web site in Application Note 12 "Avoiding Evaporation".

## 8.3 Condensation Inside the Stage Top Incubator

Check the temperature of the chamber (Heated Lid and Heated Plate). Make sure the humidity sensor is not in contact with the Heated Lid. In case of condensation, decrease the humidity and air-dry the incubator if necessary.

Please contact ibidi at techsupport@ibidi.com for further troubleshooting help.

<sup>&</sup>lt;sup>1</sup>If the experiment needs to be started immediately, either after placing the sample on the microscope or after closing the lid, we recommend checking the focus for 20 minutes. In the first few minutes after starting the experiment, temperature equilibration might influence the focus/z-position of the cells.





#### ibidi GmbH

Lochhamer Schlag 11 82166 Gräfelfing Germany

Toll free within Germany: Phone: 0800 / 00 11 11 28 Fax: 0800 / 00 11 11 29

International calls: Phone: +49 89 / 520 46 17 - 0 Fax: +49 89 / 520 46 17 - 59

E-Mail: info@ibidi.com www.ibidi.com



