Instructions





The micro–Insert 4 Well is a biocompatible, silicone insert ideal for long term microscopy and stem cell cultivation. The micro–Insert has a special sticky bottom that prevents leaking on smooth, dry surfaces, such as microscopy slides and dishes. Adherent or non–adherent cells placed into the small wells sink to the bottom and can be analyzed microscopically. The wells provide rectangular digital format (4:3) for covering the entire area by CCD cameras on microscopes. Tapered walls provide excellent optical access to cells located near the edges of the well. After cell attachment, the micro-Insert can be removed by using sterile tweezers without leaving a residue.

Material

The micro–Insert 4 Well is manufactured from biocompatible silicone. Although, the material is autoclavable and compatible to alcohols, reuse is not recommend.

Please note! When using an ibidi μ -Dish, μ -Slide or μ -Plate, make sure that the ibidi Polymer Coverslip is compatible with the immersion oil you intend to use. See page 2 for the list of compatible oils.

Geometry

Geometry of the micro–Insert 4 Well			
Number of wells	4		
Well dimensions	$2.0 \text{ mm} \times 1.5 \text{ mm}$		
Growth area per well	0.03 cm^2		
Coating area per well	0.23 cm^2		
Insert diameter	12 mm		
Height	4.2 mm		
Volume small well	10 µl		
Volume complete insert	150 µl		

We recommend using the micro–Insert 4 Well in ibidi μ –Dishes, μ –Slide 2 Well, 6–well plates, 12–well plates or large Petri dishes. They can also be used on sterile glass coverslips or glass slides.

Shipping and Storage

The μ -Slides, μ -Dishes and μ -Plates are sterilized and welded in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is listed in the following table.

Conditions		
Shipping conditions	Ambient	
Storage conditions	RT (15-25°C)	

Shelf Life of Different Surfaces			
ibiTreat, Glass Bottom, ESS	36 months		
Collagen, Poly-L-Lysine	18 months		

Surfaces and Coatings

The micro–Inserts can be transferred to any flat, clean, and dry surface. Use sterile tweezers for transfer and gently push the Insert in place. Note that only the bottom side is sticky. The micro–Insert 4 Well does not work on wet or damp surfaces. Uneven or dusty surfaces might also be problematic.

Tip:

Check whether the micro–Insert is completely attached to the surface. If not, gently push the Insert.


Seeding Cells

- Prepare cell suspension. Depending on your cell type application of a $1 7 \times 10^5$ cells/ml should result in a confluent cell layer within 24 hours.
- Add 10 µl of the cell suspension into each well. Avoid shaking as this will result in an inhomogeneous cell distribution.
- Incubate at 37°C and 5 % CO₂.
- If desired, fill the outer area with cell suspension or cell medium.
- Conduct your experiment.





micro-Insert Selection Guide



Instructions



Ordering Information

The micro–Insert is available with different well geometries and in various product versions. micro-Insert in μ –Dish^{35mm, high}

Cat. No.	Description
80406	micro–Insert 4 Well in μ –Dish ^{35mm, high} ibiTreat: ready to use, tissue culture treated, sterilized
80486	micro–Insert 4 Well FulTrac in μ–Dish ^{35mm, high} ibiTreat : ready to use, tissue culture treated, sterilized

25 micro–Inserts for self insertion



Cat. No.	Description
80409	25 micro–Inserts 4 Well for self insertion: in a 10 cm transport dish sterilized
80489	25 micro–Inserts 4 Well FulTrac for self insertion: in a 10 cm transport dish sterilized



For research use only!

Further technical specifications can be found at www.ibidi.com. For questions and suggestions please contact us by e-mail *info@ibidi.de* or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany. © ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.