



Steelco Spa – an Italian company belonging to the German Miele Group, offers a wide range of washer and autoclave models, allowing for the implementation of very complex washroom and sterilization projects. Along with the equipment, Animalab provides technical inspections as well as warranty and post-warranty service.



WASHERS FOR LABORATORY

- Equipped with various injection nozzle systems for effective cleaning of test tubes, pipettes, flasks, cylinders etc.
- Compact models with various capacities: underbench, passthrough, also with glass doors
- Various machine configurations
- Optionally equipped with a water softener, additional detergent pumps and a disinfection function
- A very large selection of inserts, baskets and trays

MODELS: LAB 500 | 600 | 610 | 640 | 660 | 680 | 900 | 1000 | 1500 | DRYING OPTION



WASHERS FOR VIVARIUM

MODELS FOR SMALLER VIVARIUM:
 AC 1200 | 1300 | 1400 | 1500

- For animal cages and bottles
- Special equipment and cycle options for washing aquatic animal tanks and rabbit cages
- Low cycle costs through reduced water and energy consumption

RACK & TUNNEL WASHERS: AC 7500 | 8500 | 3500

- High loading volume with short cycle times
- As standard, two detergent dispensers enabling setting various washing programs
- Can be used as a decontamination chamber through the option of a chemical fogging generator
- Automatic doors



ASK ABOUT SYSTEMS FOR AUTOMATIC BOTTLE EMPTYING AND FILLING AS WELL AS REMOVING AND PLACING THE BOTTLE CAPS

Why choose Steelco Miele Group washers?

- Passthrough options for higher biosafety standards BSL2, BSL3, SPF
- Very wide possibilities of washer configurations with many components and additional equipment
- Rinsing circuit separate from washing circuit in every Steelco washer model
- All internal chambers made of high quality AISI 316 steel as standard
- Temperature measurement probe PT1000 as standard

ASK ABOUT THE SERVICE OF WASHERS: INSPECTIONS, WARRANTY AND POST-WARRANTY REPAIRS

Autoclaves for scientific, industrial, pharmaceutical laboratories and animal facilities

VERTICAL AUTOCLAVES

- Compact autoclaves
- A cost effective solution for laboratory use

VS LD MODELS WITH CAPACITIES OF 90 L AND 125 L



MEDIUM AND HIGH CAPACITY AUTOCLAVES

- Pass-through options available
- The door opens manually or automatically
- AISI 316L stainless steel chamber with high efficiency heat transfer full jacket system for even heat distribution
- Doors internally made of stainless steel AISI 316L with no welding
- Heat insulation with durable reusable cover for fast service maintenance at a lower cost
- Load control probe included as standard
- Colour touch screen HMI with multi level user access system
- Built-in printer for process recording
- Non proprietary widely available durable spare parts
- ECO cool down water saving system included as standard
- Designed for BSL2, BSL3 and BSL4 laboratories



VS L MODELS WITH CAPACITIES FROM 250 L TO 920 L | BULK MODELS WITH CAPACITIES FROM 905 L TO 9 070 L

Take advantage of the Animalab experience!

- Animalab's professional consultancy in the field of project preparation, connection and installation of Steelco systems
- Full warranty and post-warranty service of washers and autoclaves
- Offer of individualized „turnkey” design solutions
- Reference list available on request