

# Aspen products from TAPVEI®

– the most natural laboratory equipment in vivarium

*Populus tremula aspen wood is a hard species, which is crucial in products for laboratory animals. It has a great influence on the parameters of moisture and ammonia absorption, as well as its practical dust-free nature. Nesting materials are extremely soft and resilient. Toys that diversify the living environment effectively activate the behavior of animals and satisfy the natural gnawing instinct, without affecting the results of experiments in the field of cytotoxicity.*


*All TAPVEI products are autoclavable and biodegradable. Their quality is ensured by ISO 9001: 2015, ISO 14001: 2015 and ISO 45001: 2018 standards. TAPVEI aspen trees are grown in areas controlled and certified by the Forest Stewardship Council (FSC).*



## Aspen bedding

*Ideal for automatic bedding dispenser. Recommended for isolators and IVC systems.*

- Moisture content ca. 10%
- Relative density 180-200 kg/m<sup>3</sup>
- Absorbing capacity 260%
- Dust content (particles < 0.25mm), less than 0.1%
- Binding of ammonium 26mg/L



## Refuges

- For mice and rats
- With solid floor recommended for wire-bottom cages
- In the shape of tunnels, houses



## Arcades

- Easier to clean - they do not have a floor
- In sizes for mice and rats



## Nesting material

- Moisture content ca. 10%
- One 3 kg bag is sufficient for up to 600 cages
- Available in 3 sizes

## Gnawing bricks

- Available in different sizes and shapes
- Reusable after washing and autoclaving

## Other accessories:

- Stairs
- Labirynths
- Cubes
- Tunnels
- Balls
- Houses



Cat. no	Name	Size (mm)	Packing
---------	------	-----------	---------

## Bedding

2HP10KG	Bedding	2x2x1	Paper bag 10 kg, Pallet (33 bags / pallet)
4HK10KG	Bedding	5x5x1	

## Nesting materials

PM90L	Nesting material for mice	ca. 2x200	Paper bag 90 L (ca. 3 kg), Pallet (21 bags / pallet)
PM90L/R	Nesting material for rats	ca. 3x200	
PM90L/2R	Nesting material for mice and rats	ca. 6x200	

## Gnawing bricks

TAPV-AB-S	Bricks S	50x10x10	1000 or 2000 pcs./pkg.
TAPV-AB-M	Bricks M	100x20x20	200 pcs./pkg.
TAPV-AB-L	Bricks L	200x43x43	50 pcs./pkg.
TAPV-AB-RABBIT	Bricks for rabbits	65x58x23	100 pcs./pkg.
TAPV-AB-BALL	Balls	Ø30	500 pcs./pkg.

## Refuges

TAPV-TUBE-S	Refuges S for mice	100x75x75	30 pcs./pkg.
TAPV-TUBE-M	Refuges M for rats	198x105x105	6 pcs./pkg.
TAPV-ARCADE14	Arcade 14 for mice	198x105x90	12 pcs./pkg.
TAPV-ARCADE17	Arcade 17 for rats	180x140x107	8 pcs./pkg.
TAPV-TUNEL	Tunnel	100x40x40, hole Ø30	100 pcs./pkg.
TAPV-LABIRYNTH	Labirynth	150x43x10, 4 holes Ø30	50 pcs./pkg.
TAPV-CUBE	Cube	60x60x60, hole Ø30	72 pcs./pkg.
TAPV-STAIRS	Stairs	75x20x208	40 pcs./pkg.
TAPV-CABIN-MOUSE	Triangular mouse cabin	155x92x85, hole Ø30	24 pcs./pkg.
TAPV-CABIN-MOUSE-STAIRS	Triangular mouse cabin with stairs	155x92x85, hole Ø30	24 pcs./pkg.
TAPV-HOUSE-MOUSE	Square house for mice	110x110x67, hole Ø40	20 pcs./pkg.
TAPV-HOUSE-GUINEA PIG	Square house for guinea pigs	220x178x138, entrance 70x90	2 pcs./pkg.
TAPV-CORNER15	Corner 15	105x75x80	60 pcs./pkg.
TAPV-CORNER15-ST	Corner 15 with stairs	105x75x80	60 pcs./pkg.