

Comparison of RoVent® systems



	RoVent® Junior	RoVent® Advanced	SomnoSuite® with RoVent®
Volume ventilation	✓	✓	✓
Pressure ventilation	✓	✓	✓
Pulse ventilation mode	✗	✓	✓
PEEP	✗	✓	✓
Assist mode	✗	✓	✓
Sigh breath	✗	✓	✓
I:E ratio	✓	✓	✓
Tidal volume — measured	✓	✓	✓
Auto set-up	Limited	✓	✓
User profiles	✗	✓	✓
Alarms (parameter)	✗	✓	✓
Trigger in	✗	✓	✓
Analog out	✗	✓	✓
Digital out	✗	✓	✓
Sync out	✗	✓	✓
ADI integration	✗	✓	✓
Data output	✗	✓	✓
Internal data logging	✗	✓	✓
Realtime clock	✗	✓	✓
Temperature monitor and homeothermic control module	✗	RightTemp® included	RightTemp® included



www.animalab.eu
info@animalab.eu