

When offering equipment for animal facilities, we use the extensive network of our suppliers and their specialization in the life science industry. We deliver ready-made workstations as well as design complex, turnkey systems to facilitate workflow management in an animal facility. Animalab also provides all equipment service for vivarium and laboratory.

### Phantom2™

- ISO Class 3 airflow
- Increased protection with Active Barrier Containment
- Ease of use with best in class ergonomic design
- More usable space with increased workspace area



Models	Working area dimensions [mm]	
	width	depth
Panthom2 Mini	724	635
Panthom2 Standard	978	635
Panthom2 Long	1232	635



### SychemCHANGE

- ISO Class 3 airflow
- H14 HEPA filters
- User-friendly touch screen interface
- Ergonomic, electrically driven legs
- Chemical and impact resistant worktop made of TRESPA®



Models	Working area dimensions [mm]	
	width	depth
SychemCHANGE level	1208	666
SychemCHANGE dipped	1208	712
SychemCHANGE XL level	1508	666
SychemCHANGE XL dipped	1508	712



- H14 HEPA filters
- Double ventilated table tops prevent any cross contamination between right and left side
- Super silent operation with noise level <52 dB(A)!


### nSAFE Mobile cage changing unit

Models	Working area dimensions [mm]	
	width	depth
nSAFE 1300	1045	677
nSAFE 1600	1345	677

### Bedding management in vivarium

Thank to cooperation with Dustcontrol we design variety of supply and disposal options in your vivarium with wide range of custom-made solutions.

- Easy supply to the dispensing silo with BigBags
- Interim storage of clean bedding in the silos
- Disposal of contaminated animal bedding into a feed hopper directly at the washer and other options
- Bedding disposal into a vacuum sealed lorry container
- Technical central unit with dust separator, vacuum producer and controls



## Bedding dispensers



### SychemFILL



### SychemFILL XL



- Unparalleled dispensing accuracy with a tolerance of  $\pm 4\%$
- Compatible with all free-flowing rodent bedding, with sizes ranging from 0,9–10 mm, as well as all commonly used wood and paper-based bedding (including corn cob).
- Custom-built, automatic bedding filling station – SychemFILL XL



### FlexFlow™

- For wide variety of bedding types
- Half the footprint of typical bedding dispensers with equal throughput capabilities
- Accurately and consistently dispenses exact amounts of bedding, up to two cages at a time

## Bedding disposal stations



### SychemWASTE

- Triple stage filtration allows quiet operating
- Double chute & sliding worksurface
- Available as SychemWASTE XL with macerator and conveyor system, to reduce bin change frequency and increases cage throughput.



### Phantom BDS

- Energy efficient – uses 62% less energy
- Increased staff protection with Container Particulate Reclamation System
- Ease of use with best in class ergonomic design
- Twice the illumination



## Air showers, material transfer hatches, pass boxes

We prepare customized projects of all pass-through solutions for vivarium. Air showers can be straight through passage, angle passage, and T passage (with 3 doors). Material transfer hatches have also many options like finishing type, HEPA filtering, DOP test connection, UV or  $H_2O_2$  sterilization etc. The overall dimensions can be adjusted for facility's requirements.

**Don't hesitate contact us to prepare installation estimate.**

