

Aquaneering Product Catalog

Table of Contents

Zebrafish Tank Configuration Options	1
Zebrafish Aquatic Housing Systems - Stand Alone with Self-Contained Filtration	2-3
Individual Light Cycle Options	4
Zebrafish Aquatic Housing Systems - Modular for Central Filtration	5
M-Rack - Build As You Grow	6-9
Large Aquaria for Aquatic Species	10-12
Zebrafish Polycarbonate Tanks	13
Breeding/Egg Collections/Hatchery Tanks	14-15
Replacement Parts and Accessories for Zebrafish Racks	16-17
Enhanced High Density Xenopus Aquatic Housing Systems	18
Standard Xenopus Aquatic Housing Systems and Accessories	19-20
Xenopus Aquatic Housing Systems - Modular for Central Filtration	21
Xenopus Frog Tanks and Accessories	22
Filters and Media	23-24
Reverse Osmosis	25-26
pH and Salinity Control	27-28
Storage Tanks and Timed Water Exchange	29-30
Monitoring and Alarms	31-32
Monitoring Accessories	33
Ultraviolet Sterilizers/Commercial UV Sterilizer/ Replacement Parts	34-37
Acrylic, Fiber Optic, and LED Lighting	38
System Chillers	39
System Heaters	40-42
Pumps	43-46
Supplies (Nets, Salt, Carbon, Test Kits, Tubing, Clamps, Bulkheads, & Suction Screens)	47-49




Prices subject to change without notice 07/2015

AQUANEERING

7960 Stromesa Court
San Diego, CA 92126 USA
1-888-689-9326 858-578-2028
info@aquaneering.com
www.aquaneering.com

AQUANEERING ZEBRAFISH TANK CONFIGURATION OPTIONS

Tank configuration is based on rack width and number of shelves.



ZEBRAFISH TANK SIZES




0.8-liter tank




1.4-liter tank




1.8-liter tank



2.8-liter tank



6.0-liter tank



9.5-liter tank

Zebrafish Aquatic Housing Systems

STANDARD ZEBRAFISH RACK & TANK SPECIFICATIONS	
	<p>Rack Specifications:</p> <ul style="list-style-type: none"> • Square T-316 Stainless Steel Construction • Removable Water Manifold, Drain Gutter, Plumbing and Connections • Water Manifolds with High Quality Flow Valves • Low Water Level Safety Shutoff
	<p>Tank Features:</p> <ul style="list-style-type: none"> • Self Cleaning Design • Autoclavable Polycarbonate • Lids and Baffles

Stand Alone Racks

Part No.#	Description	Price
ZTEDU235	<p>Zebrafish Housing, Single Sided, Two Shelf Stand Alone Education Rack 35.0" Wide x 24.0" Deep x 61" Tall 2 GPM Self-Contained Filtration System:</p> <ul style="list-style-type: none"> • Fluidized Bed Biological Filter • Carbon Finish Filter • 8-Watt UV Sterilizer • 800 Watt Titanium Electronically-Controlled Heater <p>Inline Carbon Prefilter Low Profile, Heavy Duty Plastic Cart with Casters Top Row of Tanks at Average Eye-Level (2) Quick Connect Fittings for Inlet Water (3/8") andp Overflow (3/4") Tanks with Lids and Baffles:</p> <ul style="list-style-type: none"> • 8 1.8-Liter Tanks • 5 2.8-Liter Tanks • 6 1.0-Liter Crossing Tanks 	Call for Pricing
ZT350	<p>Zebrafish Housing, Single Sided, Three Shelf Stand Alone Rack 50.5" Wide x 16.0" Deep x 52" Tall (58" Tall with UV Sterilizer) 3 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing
ZS560	<p>Zebrafish Housing, Single Sided, Five Shelf Stand Alone Rack 60.0" Wide x 14.0" Deep x 82" Tall (88" Tall with UV Sterilizer) 6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 40-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing
ZD560	<p>Zebrafish Housing, Double Sided, Five Shelf Stand Alone Rack 60.0" Wide x 26.0" Deep x 82" Tall (88" Tall with UV Sterilizer) 12 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 80-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing



Tank Configuration & Rack Height

- Tank Configuration to be Selected by Customer.
- Rack Height Will Vary Based on Shelf Spacing.

Stand Alone Racks Continued



Part No.#	Description	Price
ZS660	Zebrafish Housing, Single Sided, Six Shelf Stand Alone Rack 60.0" Wide x 14.0" Deep x 92.0" Tall (M3 89" Tall) (98" Tall with UV Sterilizer) 6 GPM Self Contained Filtration System: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 40-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing
ZD660	Zebrafish Housing, Double Sided, Six Shelf Stand Alone Rack 60.0" Wide x 26.0" Deep x 92" Tall (M3 89" Tall) (98" Tall with UV Sterilizer) 12 GPM Self Contained Filtration System: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 80-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing
ZS760	Zebrafish Housing, Single Sided, Seven Shelf Stand Alone Rack 60.0" Wide x 14.0" Deep x 89" Tall (95" Tall with UV Sterilizer) 6 GPM Self Contained Filtration System: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 40-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing
ZD760	Zebrafish Housing, Double Sided, Seven Shelf Stand Alone Rack 60.0" Wide x 26.0" Deep x 94" Tall 12 GPM Self Contained Filtration System: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 80-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater 	Call for Pricing

Tank Configuration & Rack Height

- Tank Configuration to be Selected by Customer.
- Rack Height Will Vary Based on Shelf Spacing.

Light Enclosures for Complete Rack



Part No.#	Description	Price
Light Enclosure Option:	All rack models are available with Aquaneering light enclosures. The roll up door of black ABS plastic retracts into the hood and out of the way, leaving unobstructed passageways for lab personnel. Fiber-optic lighting keeps electricity outside of the enclosure for both lab personnel safety and reliable temperature control without the need for additional water cooling equipment.	Call for Pricing
ZLHT350	Roll up Light Enclosure for Table Top Three Shelf	Call for Pricing
ZLHS560	Roll up Light Enclosure for Single Five Shelf	Call for Pricing
ZLHS660	Roll up Light Enclosure for Single Six Shelf	Call for Pricing

Light Cover for Education Rack

Part No.#	Description	Price
ZLHCCEDU	Light-tight Cover for Education Rack Black Canvas Cover Front with Full Length Velcro Closure on Each Side for Easy Access	Call for Pricing

Light Boxes for Individual Shelves




Part No.#	Description	Price
Light Box Option:	All rack models are available with Aquaneering Light Box Inserts. Feature LED Lighting, Tailgate Style Doors and Timer for Automatic Control of Breeding Light Cycle.	Call for Pricing
ZLB1060	Light Box Insert for 1.4, 1.8, 2.8, and 6.0L Rack Shelves with Tailgate Style Door and Independent Light Cycle for 60" Wide Shelf	Call for Pricing
ZLB1260	Light Box Insert for 9.5L Rack Shelves with Tailgate Style Door and Independent Light Cycle for 60" Wide Shelf	Call for Pricing

Note: For 230V 50Hz European Plug, add "EU" to end of part number.

Epoxy Powder Coat Finish

Part No.#	Description	Price
RHPCWHT-S	Rack Hardware, Powder Coat Finish, White, Single Sided Rack (Up to 75"W, 20"D, 94"H)	Call for Pricing
RHPCWHT-D	Rack Hardware, Powder Coat Finish, White, Double Sided Rack (Up to 75"W, 30"D, 94"H)	Call for Pricing

Modular Systems Requiring Central Filtration



Part No.#	Description	Price
ZAS Series		
Single sided racks for use with Aquaneering central filtration system		
ZAS550	Zebrafish Housing, Single Sided, Five Shelf Rack for Central Filtration 50.5" Wide x 14.0" Deep x 82" Tall	Call for Pricing
ZAS650	Zebrafish Housing, Single Sided, Six Shelf Rack for Central Filtration 50.5" Wide x 14.0" Deep x 92" Tall (M3 Tanks 89" T)	Call for Pricing
ZAS560	Zebrafish Housing, Single Sided, Five Shelf Rack for Central Filtration 60.0" Wide x 14.0" Deep x 82" Tall	Call for Pricing
ZAS660	Zebrafish Housing, Single Sided, Six Shelf Rack for Central Filtration 60.0" Wide x 14.0" Deep x 92" Tall	Call for Pricing
ZAS760	Zebrafish Housing, Single Sided, Seven Shelf Rack for Central Filtration 60.0" Wide x 14.0" Deep x 94" Tall (2 Rows 0.8L Tanks)	Call for Pricing



Part No.#	Description	Price
ZAD Series		
Double sided racks for use with Aquaneering central filtration system		
ZAD550	Zebrafish Housing, Double Sided, Five Shelf Rack for Central Filtration 50.5" Wide x 26.25" Deep x 82" Tall	Call for Pricing
ZAD650	Zebrafish Housing, Double Sided, Six Shelf Rack for Central Filtration 50.5" Wide x 26.25" Deep x 92" Tall (M3 Tanks 89" T)	Call for Pricing
ZAD560	Zebrafish Housing, Double Sided, Five Shelf Rack for Central Filtration 60.0" Wide x 26.25" Deep x 82" Tall	Call for Pricing
ZAD660	Zebrafish Housing, Double Sided, Six Shelf for Central Filtration 60.0" Wide x 26.25" Deep x 92" Tall (M3 Tanks 89" T)	Call for Pricing
ZAD760	Zebrafish Housing, Double Sided, Seven Shelf Rack for Central Filtration 60.0" Wide x 26.25" Deep x 94" Tall (2 Rows 0.8L Tanks)	Call for Pricing

Tank Configuration & Rack Height

- Tank Configuration to be Selected by Customer.
- Rack Height Will Vary Based on Shelf Spacing.

M-Rack Modular - Build As You Grow

STANDARD ZEBRAFISH RACK & TANK SPECIFICATIONS	
	<p>Rack Specifications:</p> <ul style="list-style-type: none"> • Square T-316 Stainless Steel Construction • Removable Water Manifold, Drain Gutter, Plumbing and Connections • Water Manifolds with High Quality Flow Valves • Low Water Level Safety Shutoff
	<p>Tank Features:</p> <ul style="list-style-type: none"> • Self Cleaning Design • Autoclavable Polycarbonate • Lids and Baffles

Build as You Grow Racks



Part No.#	Description	Price
ZTFS003-50-115	<p>"Build As You Grow" Zebrafish Housing, Base Unit For Up to 3 Shelves 50.5" (1283mm) Wide x 16" (406mm) Deep x 17" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 3 GPM (11 LPM) Self-Contained Filtration System with Lid:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer • Stainless Steel Electronically-Controlled Heater <p>For 110/115/120 Voltage, 60 Hz</p>	Call for Pricing
ZTFS003-50-230	<p>"Build As You Grow" Zebrafish Housing, Base Unit For Up to 3 Shelves 50.5" (1283mm) Wide x 16" (406mm) Deep x 17" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 3 GPM (11 LPM) Self-Contained Filtration System with Lid:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer • Stainless Steel Electronically-Controlled Heater <p>For 220/230/240 Voltage, 50 Hz</p>	Call for Pricing
ZTFS003-60-115	<p>"Build As You Grow" Zebrafish Housing, Single Shelf Base Unit 60" (1283mm) Wide x 16" (406mm) Deep x 17" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 3 GPM (11 LPM) Self-Contained Filtration System with Lid:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer • Stainless Steel Electronically-Controlled Heater <p>For 110/115/120 Voltage, 60 Hz</p>	Call for Pricing
ZTFS003-60-230	<p>"Build As You Grow" Zebrafish Housing, Single Shelf Base Unit 60" (1283mm) Wide x 16" (406mm) Deep x 17" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 3 GPM (11 LPM) Self-Contained Filtration System with Lid:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer • Stainless Steel Electronically-Controlled Heater <p>For 220/230/240 Voltage, 50 Hz</p>	Call for Pricing

Tank Configuration & Rack Height
<ul style="list-style-type: none"> • Tank Configuration to be Selected by Customer. • Rack Height Will Vary Based on Shelf Spacing.

Build As You Grow Racks Continued

Part No.#	Description	Price
ZTFS006-50-115	"Build As You Grow" Zebrafish Housing, Base Unit For Up to 5 Shelves 50.5" (1283mm) Wide x 16" (406mm) Deep x 22" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 6 GPM (22 LPM) Self-Contained Filtration System with Lid: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Finish Filter • 40-Watt UV Sterilizer • Stainless Steel Electronically-Controlled Heater For 110/115/120 Voltage, 60 Hz	Call for Pricing
ZTFS006-50-230	"Build As You Grow" Zebrafish Housing, Base Unit For Up to 5 Shelves 50.5" (1283mm) Wide x 16" (406mm) Deep x 22" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 6 GPM (22 LPM) Self-Contained Filtration System with Lid: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 40-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater • (1) Quiet One 4000 Pump 	Call for Pricing
ZTFS006-60-115	"Build As You Grow" Zebrafish Housing, Base Unit For Up to 5 Shelves 60" (1524mm) Wide x 16" (406mm) Deep x 22" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 6 GPM (22 LPM) Self-Contained Filtration System with Lid: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Finish Filter • 40-Watt UV Sterilizer • Stainless Steel Electronically-Controlled Heater For 110/115/120 Voltage, 60 Hz	Call for Pricing
ZTFS006-60-230	"Build As You Grow" Zebrafish Housing, Base Unit For Up to 5 Shelves 60" (1524mm) Wide x 16" (406mm) Deep x 22" (432mm) Tall Square T-316 Stainless Steel Powder Coated Construction Low Water Level Safety Shutoff 6 GPM (22 LPM) Self-Contained Filtration System with Lid: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Dual Carbon Filters • 40-Watt UV Sterilizer • Stainless Steel Electronically Controlled Heater • (1) Quiet One 4000 Pump • (1) Quiet One 6000 Pump For 220/230/240 Voltage, 50 Hz	Call for Pricing

Tank Configuration & Rack Height

- Tank Configuration to be Selected by Customer.
- Rack Height Will Vary Based on Shelf Spacing.

Build As You Grow Shelving Options

Part No.#	Description	Price
ZAT150-080	<p>"Build As You Grow" Add-on Shelf for 0.8L Tanks 50.5" (1283mm) Wide x 16" (406mm) Deep x 8" (203mm) Tall Square T-316 Stainless Steel Powder Coated Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with 15 High Quality Flow Valves Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT160-080	<p>"Build As You Grow" Add-on Shelf for 0.8L Tanks 60" (1283mm) Wide x 16" (406mm) Deep x 8" (203mm) Tall Square T-316 Stainless Steel Powder Coated Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with 18 High Quality Flow Valves Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT150-M3	<p>"Build As You Grow" Add-on Shelf for 1.4, 1.8, 2.8, & 6.0L Tanks 50.5" (1283mm) Wide x 16" (406mm) Deep x 11" (279mm) Tall Square T-316 Stainless Steel Powder Coated Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with 15 High Quality Flow Valves Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT160-M3	<p>"Build As You Grow" Add-on Shelf for 1.4, 1.8, 2.8, & 6.0L Tanks 60" (1283mm) Wide x 16" (406mm) Deep x 11" (279mm) Tall Square T-316 Stainless Steel Powder Coated Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with 18 High Quality Flow Valves Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT150-950	<p>"Build As You Grow" Add-on Shelf for 9.5L Tanks 50.5" (1283mm) Wide x 16" (406mm) Deep x 12" (305mm) Tall Square T-316 Stainless Steel Powder Coated Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with 12 High Quality Flow Valves Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT160-950	<p>"Build As You Grow" Add-on Shelf for 9.5L Tanks 60" (1283mm) Wide x 16" (406mm) Deep x 12" (305mm) Tall Square T-316 Stainless Steel Powder Coated Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with 12 High Quality Flow Valves Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing

Build As You Grow Shelving with Light Box Options

Part No.#	Description	Price
ZAT150-M3LB10	<p>"Build As You Grow" Add-on Shelf with Light Box for M3 Tanks 50.5" (1283mm) Wide x 16" (406mm) Deep x 11" (279mm) Tall Square T-316 Stainless Steel Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with High Quality Flow Valves Light Box Insert for 1.4, 1.8, 2.8 or 6.0L tanks with Tailgate Style Door and Independent Light Cycle Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT160-M3LB10	<p>"Build As You Grow" Add-on Shelf with Light Box for M3 Tanks 60" (1283mm) Wide x 16" (406mm) Deep x 11" (279mm) Tall Square T-316 Stainless Steel Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with High Quality Flow Valves Light Box Insert for 1.4, 1.8, 2.8 or 6.0L tanks with Tailgate Style Door and Independent Light Cycle Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT150-950LB12	<p>"Build As You Grow" Add-on Shelf with Light Box for 9.5L Tanks 50.5" (1283mm) Wide x 16" (406mm) Deep x 12" (305mm) Tall Square T-316 Stainless Steel Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with High Quality Flow Valves Light Box Insert for 9.5L tanks with Tailgate Style Door and Independent Light Cycle Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing
ZAT160-950LB12	<p>"Build As You Grow" Add-on Shelf with Light Box for 9.5L Tanks 60" (1283mm) Wide x 16" (406mm) Deep x 12" (305mm) Tall Square T-316 Stainless Steel Construction Removable Water Manifold, Drain Gutter, Plumbing and Connections Water Manifolds with High Quality Flow Valves Light Box Insert for 9.5L tanks with Tailgate Style Door and Independent Light Cycle Tank Features: Self-Cleaning Design Autoclavable Polycarbonate Tanks with Lids and Baffles (included with rack price):</p>	Call for Pricing

Large Aquaria Housing for FG Series

Standard Aquatic Rack Specifications	
Rack Specifications: <ul style="list-style-type: none"> • Square T-316 Stainless Steel Construction • Removable Water Manifold, Drain Gutter, Plumbing and Connections • Water Manifolds with High Quality Flow Valves • Low Water Level Safety Shutoff 	

Stand Alone Racks

Part No.#	Description	Price
FG Series Stand Alone Rack		
Racks for Glass Aquaria and 6 GPM Self Contained Recirculation System (Unless Otherwise Stated). Tanks Sold Separately.		
FG350	Large Aquaria Housing, Single Sided, Three Shelf Stand Alone Rack 50.5" Wide x 22.5" Deep x 82.0" Tall 3 GPM Self Contained Filtration System: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer 	Call for Pricing
FG360	Large Aquaria Housing, Single Sided, Three Shelf Stand Alone Rack 60.0" Wide x 22.5" Deep x 82.0" Tall	Call for Pricing
FG375	Large Aquaria Housing, Single Sided, Three Shelf Stand Alone Rack 75.0" Wide x 22.5" Deep x 82.0" Tall	Call for Pricing
Options include Air lines and Shelf Mounted Electrical Outlets with In-use Covers.		

Part No.#	Description	Price
FP Series Stand Alone Rack		
Racks for Glass Aquaria and 6 GPM Self Contained Recirculation System (Unless Otherwise Stated). Tanks Sold Separately.		
FP350	Large Aquaria Housing, Single Sided, Three Shelf Stand Alone Rack 50.5" Wide x 24.0" Deep x 82.0" Tall 3 GPM Self Contained Filtration System: <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Carbon Finish Filter • 25-Watt UV Sterilizer 	Call for Pricing
FP360	Large Aquaria Housing, Single Sided, Three Shelf Stand Alone Rack 60.0" Wide x 24.0" Deep x 82.0" Tall	Call for Pricing
FP375	Large Aquaria Housing, Single Sided, Three Shelf Stand Alone Rack 75.0" Wide x 24.0" Deep x 82.0" Tall	Call for Pricing
Options include Air lines and Shelf Mounted Electrical Outlets with In-use Covers.		



Tank Configuration & Rack Height

- Tank Configuration to be Selected by Customer.
- Rack Height Will Vary Based on Shelf Spacing.

Modular Systems Requiring Central Filtration



Part No.#	Description	Price
FAG Series Stand Alone Rack		
Racks for Glass Aquaria and 6 GPM Self Contained Recirculation System (Unless Otherwise Stated). Tanks Sold Separately.		
FAG350	Large Aquaria Housing, Single Sided, Three Shelf Rack for Central Filtration 50.5" Wide x 22.5" Deep x 82.0" Tall	Call for Pricing
FAG360	Large Aquaria Housing, Single Sided, Three Shelf Rack for Central Filtration 60.0" Wide x 22.5" Deep x 82.0" Tall	Call for Pricing
FAG375	Large Aquaria Housing, Single Sided, Three Shelf Rack for Central Filtration 75.0" Wide x 22.5" Deep x 82.0" Tall	Call for Pricing
Racks for Glass Aquaria and 6 GPM Self Contained Recirculation System (Unless Otherwise Stated). Tanks Sold Separately.		



Part No.#	Description	Price
FAP Series Stand Alone Rack		
Racks for Glass Aquaria and 6 GPM Self Contained Recirculation System (Unless Otherwise Stated). Tanks Sold Separately.		
FAP350	Large Aquaria Housing, Single Sided, Three Shelf Rack for Central Filtration 50.5" Wide x 24.0" Deep x 82.0" Tall	Call for Pricing
FAP360	Large Aquaria Housing, Single Sided, Three Shelf Rack for Central Filtration 60.0" Wide x 24.0" Deep x 82.0" Tall	Call for Pricing
FAP375	Large Aquaria Housing, Single Sided, Three Shelf Rack for Central Filtration 75.0" Wide x 24.0" Deep x 82.0" Tall	Call for Pricing
Racks for Glass Aquaria and 6 GPM Self Contained Recirculation System (Unless Otherwise Stated). Tanks Sold Separately.		

Glass Tanks for FG Series Racks

Part No.#	Description	Price
FTG02.5-6128	2.5 Gallon Glass Tank 6" x 12" x 8" Tall (Drilled)	
FTG05-5.52012	5 Gallon Glass Tank 5.5" x 20" x 12" Tall (Drilled)	
FTG10-112012	10 Gallon Glass Tank 11" x 20" x 12" Tall (Drilled)	
FTG20-222012	20 Gallon Glass Tank 22" x 20" x 12" Tall (Drilled)	
FTG30-332012	30 Gallon Glass Tank 33" x 20" x 12" Tall (Drilled)	
FTG40-442012	40 Gallon Glass Tank 44" x 20" x 12" Tall (Drilled)	
FTG50-502012	50 Gallon Glass Tank 50" x 20" x 12" Tall (Drilled)	
Other standard Glass Aquaria sizes available.		


Polycarbonate Tanks

Part No.#	Description	Price
FTP08	8 Gallon Polycarbonate Tank 11" x 21" x 12" Tall Drilled with Drain Assembly	
FTP20	20 Gallon Polycarbonate Tank 24" x 21" x 12" Tall Drilled with Drain Assembly	
FTP30	30 Gallon Polycarbonate Tank 35" x 21" x 12" Tall Drilled with Drain Assembly	


Polyethylene Tanks

Part No.#	Description	Price
FTP016	16 Liter Molded Polycarbonate Tank with Screened Rear Drain, Sliding Lid with Air and Water Openings Enclosed by Grommet to Prevent Cross Contamination from Above Tanks	
FTP050	50 Liter Molded Polyethylene Tank with Center Double Standpipe and Bottom Drain, 6" Radius Rounded Corners. Translucent Sliding Lid with Air and Water Openings Enclosed by Grommet to Prevent Cross Contamination from Above Tanks. Plugs for Exposed Holes After Tank is Removed Included.	
FTP100	100 Liter Molded Polyethylene Tank with Center Double Standpipe and Bottom Drain, 6" Radius Rounded Corners. Translucent Sliding Lid with Air and Water Openings Enclosed by Grommet to Prevent Cross Contamination from Above Tanks. Plugs for Exposed Holes After Tank is Removed Included.	
FTP800	800 Liter Molded Polyethylene Tank with Center Double Standpipe, Stand, and Bottom Drain, 6" Radius Rounded Corners. Molded Lid (same color as tank) with Air and Water Openings.	


Zebrafish Polycarbonate Tanks




Part No.#	Description	Price
ZT080T	0.8 Liter Clear Tank	
ZT080L	0.8 Liter Green Lid	
ZT080B	0.8 Liter Green Baffle	
ZT080	Complete 0.8 Liter Clear Tank Set includes: Tank, Green Lid and Green Baffle	
ZT080S400	0.8 Liter Screen - 400 Micron (optional)	
ZT080S850	0.8 Liter Screen - 850 Micron (optional)	




Part No.#	Description	Price
ZT140T	1.4 Liter Clear Tank	
ZT140L	1.4 Liter Green Lid	
ZT140B	1.4 Liter Green Baffle	
ZT140	Complete 1.4 Liter Clear Tank Set includes: Tank, Green Lid and Green Baffle	
ZT140S400	1.4 Liter Screen - 400 Micron (optional)	
ZT140S850	1.4 Liter Screen - 850 Micron (optional)	




Part No.#	Description	Price
ZT180T	1.8 Liter Clear Tank	
ZT180L	1.8 Liter Green Lid	
ZT180B	1.8 Liter Green Baffle	
ZT180	Complete 1.8 Liter Clear Tank Set includes: Tank, Green Lid and Green Baffle	
ZT180S400	1.8 Liter Screen - 400 Micron (optional)	
ZT180S850	1.8 Liter Screen - 850 Micron (optional)	



Part No.#	Description	Price
ZT280T	2.8 Liter Clear Tank	
ZT280L	2.8 Liter Green Lid	
ZT280B	2.8 Liter Green Baffle	
ZT280	Complete 2.8 Liter Clear Tank Set includes: Tank, Green Lid and Green Baffle	
ZT280S400	2.8 Liter Screen - 400 Micron (optional)	
ZT280S850	2.8 Liter Screen - 850 Micron (optional)	



Part No.#	Description	Price
ZT600T	6.0 Liter Clear Tank	
ZT600L	6.0 Liter Green Lid	
ZT600B	6.0 Liter Green Baffle	
ZT600	Complete 6.0 Liter Clear Tank Set includes: Tank, Green Lid and Green Baffle	
ZT600S400	6.0 Liter Screen - 400 Micron (optional)	
ZT600S850	6.0 Liter Screen - 850 Micron (optional)	



Part No.#	Description	Price
ZT950T	9.5 Liter Clear Tank	
ZT950L	9.5 Liter Green Lid	
ZT950B	9.5 Liter Green Baffle	
ZT950	Complete 9.5 Liter Clear Tank Set includes: Tank, Green Lid and Green Baffle	
ZT950S400	9.5 Liter Screen - 400 Micron (optional)	
ZT950S850	9.5 Liter Screen - 850 Micron (optional)	

Part No.#	Description	Price
ZTP-3000	30 Liter Acrylic Tank for Zebrafish Rack 27" x 11" x 8" Tall with Lid	
ZTP-5000	50 Liter Acrylic Tank for Zebrafish Rack 48" x 11" x 8" Tall with Lid	
ZTP-6000	60 Liter Acrylic Tank for Zebrafish Rack 55" x 11" x 8" Tall with Lid	

DC 96 Holding System

Part No.#	Description	Price
RHDC24	24-Compartment housing tank with removable dividers	
RHDC24SD	Sliding drawer unit for 24-compartment housing tank. Mounts to 0.8-liter tank shelf	
RHDC24BA	Bracket assembly for 96-compartment housing tank - Aquaneering Rack	
RHDC24SET	24-Compartment housing tank with sliding drawer unit and removable dividers.	
RHDC24PK	Plumbing Kit for one shelf of three tanks (with sliding drawers or brackets), includes bar, (3) mini tank manifolds, and (1) supply manifold.	

Breeding/Egg Collections/Hatchery Tank

Part No.#	Description	Price
ZHBS02	2.0L Brine Shrimp Hatcher with Stand	
ZHBS02H	2.0L Brine Shrimp Hatcher without Stand	
ZHBS02S	2.0L Brine Shrimp Hatcher Stand only	
ZHBS20	20L Brine Shrimp Hatcher with Stand	
ZHBS20H	20L Brine Shrimp Hatcher without Stand	
ZHBS25	25L White Polyethylene Brine Shrimp Hatcher with Stand	
ZHBS25H	25L White Polyethylene Brine Shrimp Hatcher without Stand	
ZHBS20S	20L Brine Shrimp Hatcher Stand only	
MH100GC	Heater 100 Watt Drop In, Glass, with Built-In Controller for 20L Brine Shrimp Hatcher	
ZHBS100	100L Brine Shrimp Hatcher with Powder Coated Stand Includes 2 Airstones and 10 ft. of Airline	
ZHBS100H	100L Brine Shrimp Hatcher Only	
ZHBSAP	Air Pump for 2L Brine Shrimp Hatcher-Single Outlet	
ZHBDAP	Air Pump for 2L Brine Shrimp Hatcher-Double Outlet	
ZHBSAP-20	Air Pump for 20L Brine Shrimp Hatcher-Single Outlet Adjustable Air Flow	
ZHBDAP-20	Air Pump for 20L Brine Shrimp Hatcher-Double Outlet	
ZHCT100	Crossing Tank (1.0L) Clear Polycarbonate with Clear Lid and Clear Insert	
ZHCT100T	Crossing Tank (1.0-liter) Clear Polycarbonate	
ZHCT100L	Crossing Tank Lid (1.0-liter) Clear Polycarbonate	
ZHCT100I	Crossing Tank Insert (1.0-liter) Clear Polycarbonate	
ZHCT100D	Crossing Tank Divider (1.0-liter) Clear Polycarbonate	
ZHCCTray-1015W	Crossing Rack Tray, 10" x 15" White (Single Sided Rack Shelving)	
ZHHCTray-1015B	Crossing Rack Tray, 10" x 15" Black (Single Sided Rack Shelving)	
ZHCCTray-926W	Crossing Rack Tray, 9" x 26" White (Double Sided Rack Shelving)	
ZHCCTray-926B	Crossing Rack Tray, 9" x 26" Black (Double Sided Rack Shelving)	
ZHCCTray-1224W	Crossing Rack Tray, 12" x 24" White (Extra Deep for Large Tanks)	
ZHCCTray-1224B	Crossing Rack Tray, 12" x 24" Black (Extra Deep for Large Tanks)	

Breeding/Egg Collections/Hatchery Tanks Cont.

Part No.#	Description	Price
ZHCC0510	Crossing Tank Cart with Casters & Five Trays to Hold 50 Crossing Tanks	
ZHCC1010	Crossing Tank Cart with Casters & Ten Trays to Hold 100 Crossing Tanks	
ZHCS1005	Crossing Tank Shelves - Holds (50) Crossing Tanks 26" Wide x 10" Deep x 82" Tall with 10" Half Size Trays Mounted to Wall (Tanks not included)	
ZHCS1005RM	Crossing Tank Shelves - Holds (50) Crossing Tanks 26" Wide x 10" Deep x 82" Tall with 10" Half Size Trays With Hardware to Mount to End of Double Rack (Tanks not included)	
ZHESC100	Zebrafish Embryo Spawning Chamber 20" Dia x 25" Deep x 48" Tall with Stand 100-Liter Spawning Chamber with Removable Egg Traps, Off-Cycle Lighting, Plastic Stand, and System Connection Kit	
ZTTEESC-2	Embryo Spawning System 72" Wide x 22" Deep x 48.5" Tall Two Spawning Chambers Mounted on Stainless Steel Stand (Includes Light Controllers) 3 GPM Filtration System Fluidized Bed Biofilter Biological Filtration Activated Carbon Polishing Filtration 25 Watt UV Sterilizer Stainless Steel Heater	
ZHNT2-280	Zebrafish Hatchery Set: 2 Hatchers with 2.8 Liter Lid	
ZHNT4-600	Zebrafish Hatchery Set: 4 Hatchers with 6.0 Liter Lid	
ZHNT4-950	Zebrafish Hatchery Set: 4 Hatchers with 9.5 Liter Lid	
ZHNTOR	Zebrafish Hatchery Set O-Ring	
ZHEC280	2.8 Liter Egg Collection Box	
ZHEC950	9.5 Liter Egg Collection Box	

Replacement Parts and Accessories for Zebrafish Racks

Manifolds and Replacement Parts

Part No.#	Description	Price
ZMR3	Riser Manifold Assembly (3 Shelf)	
ZMR5	Riser Manifold Assembly (5 Shelf)	
ZMR6	Riser Manifold Assembly (6 Shelf)	
ZPR007-62	Riser Pipe	
ZMV1050NV	Valve Manifold Water (10 Valve) for 50.5" Rack (Valves Not Included)	
ZMV1550NV	Valve Manifold Water (15 Valve) for 50.5" Rack (Valves Not Included)	
ZMV1650NV	Valve Manifold Water (16 Valve) for 50.5" Rack (Valves Not Included)	
ZMV0660NV	Valve Manifold Water (6 Valve) for 60" Rack (Valves Not Included)	
ZMV1260NV	Valve Manifold Water (12 Valve) for 60" Rack (Valves Not Included)	
ZMV1860NV-S	Valve Manifold Water (18 Valve) for 60" Rack - Short (Valves Not Included)	
ZMV1860NV	Valve Manifold Water (18 Valve) for 60" Rack (Valves Not Included)	
ZMV2060NV	Valve Manifold Water (20 Valve) for 60" Rack (Valves Not Included)	
ZMV1469NV	Valve Manifold Water (14 Valve) for 69" Rack (Valves Not Included)	
ZMV2160NV	Valve Manifold Water (21 Valve) for 69" Rack (Valves Not Included)	
ZMV2477NV	Valve Manifold Water (24 Valve) for 77" Rack (Valves Not Included)	
ZMV1050	Valve Manifold Water (10 Valve) for 50.5" Rack (Valves Included)	
ZMV1550	Valve Manifold Water (15 Valve) for 50.5" Rack (Valves Included)	
ZMV1650	Valve Manifold Water (16 Valve) for 50.5" Rack (Valves Included)	
ZMV0660	Valve Manifold Water (6 Valve) for 60" Rack (Valves Included)	
ZMV1260	Valve Manifold Water (12 Valve) for 60" Rack (Valves Included)	
ZMV1860-S	Valve Manifold Water (18 Valve) for 60" Rack - Short (Valves Included)	
ZMV1860	Valve Manifold Water (18 Valve) for 60" Rack (Valves Included)	
ZMV2060	Valve Manifold Water (20 Valve) for 60" Rack (Valves Included)	
ZMV1469	Valve Manifold Water (14 Valve) for 69" Rack (Valves Included)	
ZMV2160	Valve Manifold Water (21 Valve) for 69" Rack (Valves Included)	
ZMV2477	Valve Manifold Water (24 Valve) for 77" Rack (Valves Included)	
PTBLU025	Tank Supply Tubing - Blue (Flexible) 1/4" OD per ft.	
PTBLU025-100	Tank Supply Tubing - Blue (Flexible) 1/4" OD 100 ft. roll	
PTBLU025-500	Tank Supply Tubing - Blue (Flexible) 1/4" OD 500 ft. roll	
ZMV	Tank Supply Valve	
ZTLG	Lid Grommets (fits 1.4, 2.8, 9.5 Liter Tanks)	



Downspouts & Gutters

Part No.#	Description	Price
ZDA3	Downspout Assembly (3 Shelf)	
ZDA5	Downspout Assembly (5 Shelf)	
ZDA6	Downspout Assembly (6 Shelf)	
ZG35	Gutter (35" Wide for eRack)	
ZG50	Gutter (50.5" Wide Rack)	
ZG60	Gutter (60.0" Wide Rack)	

Miscellaneous

Part No.#	Description	Price
SHDSB	Downspout Scrub Brush with 72" Extension	
ZFV002	Float Valve 1/4"	
ZFV007	Float Valve 3/4"	
ZSB50	Stainless Steel 9.5 Liter Tank Support Bar for 50.5" Rack	
ZSB60	Stainless Steel 9.5 Liter Tank Support Bar for 60" Rack	
TOOL-VRJG	John Guest Valve Removal Tool	


Enhanced High Density Xenopus Aquatic Housing Systems

STANDARD Xenopus RACK & TANK SPECIFICATIONS



Rack Specifications:

- Square T-316 Stainless Steel Construction
- Removable Water Manifold, Drain Gutter, Plumbing and Connections
- Water Manifolds with High Quality Flow Valves
- Low Water Level Safety Shutoff


Stand Alone Racks



Part No.#	Description	Price
XS250-421-250	<p>Xenopus Housing, Two Shelf Stand Alone Rack 50.5" Wide x 20.0" Deep x 58" Tall (64" Tall with UV Sterilizer)</p> <p>3 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 25-Watt UV Sterilizer <p>Roto-Molded Sump Tank with Cover Slide-Out Tanks for Easy Removal</p> <p>Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 4 21-Liter Clear PETG Tanks • 2 50-Liter Clear PETG Tanks 	Call for Pricing
XS360-521-650	<p>Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 82" Tall (88" Tall with UV Sterilizer)</p> <ul style="list-style-type: none"> • Diverter Valve to connect Top Shelf to "Flow-Through" plumbing for both drain and supply, to allow connection to "Timed" water exchange. • Valves to be Labeled: Flow Thru/Recirculation <p>Slide-Out Tanks for Easy Removal</p> <p>6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 40-Watt UV Sterilizer <p>Roto-Molded Sump Tank with Cover</p> <p>Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 5 21-Liter Clear PETG Tanks • 6 50-Liter Clear PETG Tanks 	Call for Pricing
XS360-1521	<p>Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 82" Tall (88" Tall with UV Sterilizer)</p> <ul style="list-style-type: none"> • Diverter Valve to connect Top Shelf to "Flow-Through" plumbing for both drain and supply, to allow connection to "Timed" water exchange. • Valves to be Labeled: Flow Thru/Recirculation <p>Slide-Out Tanks for Easy Removal</p> <p>6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 40-Watt UV Sterilizer <p>Roto-Molded Sump Tank with Cover</p> <p>Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 15 21-Liter Clear PETG Tanks 	Call for Pricing




Stand Alone Racks Cont.





Part No.#	Description	Price
XS360-950	<p>Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 82" Tall (88" Tall with UV Sterilizer)</p> <ul style="list-style-type: none"> • Diverter Valve to connect Top Shelf to "Flow-Through" plumbing for both drain and supply, to allow connection to "Timed" water exchange. • Valves to be Labeled: Flow Thru/Recirculation <p>Slide-Out Tanks for Easy Removal 6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 40-Watt UV Sterilizer <p>Roto-Molded Sump Tank with Cover Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 9 50-Liter Clear PETG Tanks 	Call for Pricing

Standard Xenopus Aquatic Housing Systems Basic, More Economical Model.


Stand Alone Racks





Part No.#	Description	Price
XT250-01	<p>Xenopus Housing, Two Shelf Stand Alone Rack 50.5" Wide x 20.0" Deep x 53" Tall (59" Tall with UV Sterilizer)</p> <p>3 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 25-Watt UV Sterilizer <p>Roto-Molded Sump Tank Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 8 16-Liter Clear Polycarbonate Tanks 	Call for Pricing
XT250-02	<p>Xenopus Housing, Two Shelf Stand Alone Rack 50.5" Wide x 20.0" Deep x 55" Tall (61" Tall with UV Sterilizer)</p> <p>3 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 25-Watt UV Sterilizer <p>Roto-Molded Sump Tank Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 4 16-Liter Clear Polycarbonate Tanks • 2 50-Liter (White Polyethylene Tanks or Clear PETG +\$400) 	Call for Pricing
XT250-03	<p>Xenopus Housing, Two Shelf Stand Alone Rack 50.5" Wide x 20.0" Deep x 57" Tall (63" Tall with UV Sterilizer)</p> <p>3 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 25-Watt UV Sterilizer <p>Roto-Molded Sump Tank Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 4 50-Liter Clear PETG Tanks 	Call for Pricing



Stand Alone Racks Cont.



Part No.#	Description	Price
XS360-01	<p>Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 76" Tall (82" Tall with UV Sterilizer)</p> <p>6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 40-Watt UV Sterilizer <p>Roto-Molded Sump Tank</p> <p>Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 15 16-Liter Clear Polycarbonate Tanks 	Call for Pricing
XS360-02	<p>Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 83" Tall (86" Tall with UV Sterilizer)</p> <p>Heavy-Duty Stainless Steel Castors</p> <p>6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 40-Watt UV Sterilizer <p>Roto-Molded Sump Tank</p> <p>Slide-Out Tanks for Easy Removal</p> <p>Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 5 16-Liter Clear Polycarbonate Tanks • 4 50-Liter (White Polyethylene Tanks or Clear PETG +\$400) 	Call for Pricing
XS360-03	<p>Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 86" Tall (92" Tall with UV Sterilizer)</p> <p>Heavy-Duty Stainless Steel Castors</p> <p>6 GPM Self Contained Filtration System:</p> <ul style="list-style-type: none"> • Dacron Pad Mechanical Filtration • Fluidized Bed Biological Filter • Finish Filter with Activated Carbon • 40-Watt UV Sterilizer <p>Roto-Molded Sump Tank</p> <p>Slide-Out Tanks for Easy Removal</p> <p>Tanks with Lids and Standpipes:</p> <ul style="list-style-type: none"> • 6 50-Liter (White Polyethylene Tanks or Clear PETG +\$600) 	Call for Pricing



Modular Systems for Central Filtration

Part No.#	Description	Price
XAS Series		
Single sided racks for use with Aquaneering central filtration system		
XAS350	Xenopus Housing, Three Shelf Rack for Central Filtration 50.5" Wide x 20" Deep x 83" Tall Roto-Molded Sump Tank Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 4 16-Liter Clear Polycarbonate Tanks • 4 50-Liter Clear PETG Tanks 	Call for Pricing
XAS360-01	Xenopus Housing, Three Shelf Rack for Central Filtration 60.0" Wide x 20.0" Deep x 83" Tall Roto-Molded Sump Tank Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 15 16-Liter Clear Polycarbonate Tanks 	Call for Pricing
XAS360-02	Xenopus Housing, Three Shelf Rack for Central Filtration 60.0" Wide x 20.0" Deep x 83" Tall Heavy-Duty Stainless Steel Castors Roto-Molded Sump Tank Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 5 16-Liter Clear Polycarbonate Tanks • 4 50-Liter (White Polyethylene Tanks or Clear PETG +\$400) 	Call for Pricing
XAS360-03	Xenopus Housing, Three Shelf Stand Alone Rack 60.0" Wide x 20.0" Deep x 83" Tall Heavy-Duty Stainless Steel Castors Roto-Molded Sump Tank Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 4 50-Liter (White Polyethylene Tanks or Clear PETG +\$600) 	Call for Pricing
XAS360-521-650	Xenopus Housing, Three Shelf Rack for Central Filtration 60.0" Wide x 20" Deep x 82" Tall Roto-Molded Sump Tank with Cover Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 5 21-Liter Clear PETG Tanks • 6 50-Liter Clear PETG Tanks 	Call for Pricing
XAS360-950	Xenopus Housing, Three Shelf Rack for Central Filtration 60.0" Wide x 20" Deep x 82" Tall Roto-Molded Sump Tank with Cover Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 5 21-Liter Clear PETG Tanks • 6 50-Liter Clear PETG Tanks 	Call for Pricing
XAS360-1521	Xenopus Housing, Three Shelf Rack for Central Filtration 60.0" Wide x 20" Deep x 82" Tall Roto-Molded Sump Tank with Cover Slide-Out Tanks for Easy Removal Tanks with Lids and Standpipes: <ul style="list-style-type: none"> • 15 21-Liter Clear PETG Tanks 	Call for Pricing

Xenopus Frog Tanks and Accessories

Tanks feature a plunger for one-step cleaning and sliding lid that does not interfere with water supply

COLOR CODING

C Clear
W White
A Amber Tinted

Enhanced 21 Liter Tanks



Part No.#	Description	Price
XT21C	Complete 21 Liter Clear *PETG Tank Set includes: Tank, Lid and Standpipe	
XT21CT	21 Liter Clear PETG Tank	
XT21CL	21 Liter Clear PETG Lid	

Enhanced 50 Liter Tanks



Part No.#	Description	Price
XT50C-M	Complete Enhanced 50 Liter Clear *PETG Tank Set includes: Tank, Lid and Standpipe	
XT50CT-M	Enhanced 50 Liter Clear PETG Tank	
XT50CL-M	Enhanced 50 Liter Clear PETG Lid	
XTSPA-M	Standpipe Assembly (for Clear Tanks)	


Standard Xenopus Frog Tanks and Accessories

Standard 16 Liter Tanks



Part No.#	Description	Price
XT16CT	16 Liter Clear Polycarbonate Tank	
XT16CL	16 Liter Clear Polycarbonate Lid	
XT16C	Complete 16 Liter Clear Polycarbonate Tank Set includes: Tank, Lid and Bulkhead	
XT16B	XFBHD 16 Liter Bulkhead	
XT16D	16 Liter Standpipe/Drain Assembly	
XT16PS	16 Liter Polysulphone Tank Set includes: Tank, Lid, Bulkhead and Drain Assembly	

Standard 50 Liter Tanks



Part No.#	Description	Price
XT50CT	50 Liter Clear PETG Tank	
XT50CL	50 Liter Clear PETG Lid	
XT50C	Complete 50 Liter Clear *PETG Tank Set includes: Tank, Lid and Standpipe	
XT50WT	50 Liter White Translucent Polyethylene Tank	
XT50WL	50 Liter White Translucent Polyethylene Lid	
XT50W	Complete 50 Liter White Translucent Polyethylene Tank Set includes: Tank, Lid and Standpipe	
XTSPA	Standpipe Assembly (for 21 and 50 Liter White or Clear Tanks)	
XT50RG	Removable Grate for Translucent Lid	

*PETG is a bisphenol free plastic.

Carts

XS360CART	Low Profile Heavy Duty Stainless Steel Cart w/Casters for 60" Rack. Adds 1.5" to overall height of rack.	
XS250CART	Low Profile Heavy Duty Stainless Steel Cart w/Casters for 50.5" Rack. Adds 1.5" to overall height of rack.	

Stand Alone Replacement Parts

Part No.#	Description	Price
MFAP612	Prefilter Pads (for Table Top Racks) 30 Count (6x12 for 12 x 48 x 15 1/2" Sump)	
MFAP912	Prefilter Pads (for Single Racks) 30 Count (9x12 for 17 1/2 x 58 x 16 1/2" Sump)	
MFAP622	Prefilter Pads (for Double Racks) 30 Count (6x22 17 1/2 x 58 x 26" Sump)	
MFCF12	Carbon Filter (12") with Sleeve and Carbon	
MFCF12S2	Carbon Filter Sleeve (12") 2 per pkg. (for Above Sump Carbon Filter)	
MFCF12SC	Carbon Filter Sleeve (12") 10 Case	

Reusable Filters

Part No.#	Description	Price
MFAR612	Non Consumable Reusable Prefilter (for Table Top Sumps) 6" x 12"	
MFAR912	Non Consumable Reusable Prefilter (for Single Sumps) 9" x 12"	
MFAR622	Non Consumable Reusable Prefilter (for Double Sumps) 6" x 22"	

Mechanical Filter Replacement Parts

Part No.#	Description	Price
EACVA24T	Hayward Pool 3 Way Valve Electrical Actuator for use with Automatic Back-wash Filtration	
PH3WVEACT	Hayward Commercial 3 Way Valve Electrical Actuator for use with propeller bead filters.	

FSI Filters

Part No.#	Description	Price
MFVB005	5 Micron FSI Filter Bag (each)	
MFVB005C	5 Micron FSI Filter Bag (50 Case)	
MFVB010	10 Micron FSI Filter Bag (each)	
MFVB010C	10 Micron FSI Filter Bag (50 Case)	
MFVB025	25 Micron FSI Filter Bag (each)	
MFVB025C	25 Micron FSI Filter Bag (50 Case)	
MFVB050	50 Micron FSI Filter Bag (each)	
MFVB050C	50 Micron FSI Filter Bag (50 Case)	
MFVB100	100 Micron FSI Filter Bag (each)	
MFVB100C	100 Micron FSI Filter Bag (50 Case)	
MFVA100	FSI Filter Vessel (X100) Complete Housing	
MFVBASK	FSI Filter Housing Inner Basket	
MFVL	Replacement Lid - No O-Ring for FSI Filter Vessel	
MFVG	Replacement O-Ring for FSI Filter Vessel	

Fluidized Bed Media for Stand Alone Racks

Part No.#	Description	Price
MFFMG10	Fluidized Bed Glass Media Washed 10 lb bag	
MFFMG16	Fluidized Bed Glass Media Washed 16 lb bag (Education Rack)	
MFFMG25	Fluidized Bed Glass Media Washed 25 lb bag (Table Top Zebrafish Rack)	
MFFMG50	Fluidized Bed Glass Media Washed 50 lb bag (Single Sided Zebrafish Rack)	

Fluidized Bed & Mechanical Backwashing Filter Media

Part No.#	Description	Price
AFMM	Advanced Filter, Mechanical Media Washed (50 lb bag)	
AFMS06	Advanced Filter Media, Silica #6 (50 lb bag)	
AFMS60	Advanced Filter Media, Silica #60 (100 lb bag)	
AFMS90	Advanced Filter Media, Silica #90 (100 lb bag)	

Water Filters, Carbon and Sediment

Part No.#	Description	Price
MFH10	Carbon Pre-Filter Housing - 10" (Req 1 - RO Carbon Block Cartridge 10" (MROTCB10-2.5-10))	
MFH10D	Carbon Pre-Filter Double Housing - 10" (Req 2 - RO Carbon Block Cartridges 10" (MROTCB10-2.5-10))	
MFH10DTP	Carbon Pre-Filter Double Housing with Test Port - 10" (Req 2 - RO Carbon Block Cartridges 10" (MROTCB10-2.5-10))	
MFH10T	Carbon Pre-Filter Triple Housing - 10" (Req 2 - RO Carbon Block Cartridges 10" (MROTCB10-2.5-10) and 1 - RO Sediment Cartridge 10" (MROTSC10-2.5-5))	
MFH10TP	Carbon Pre-Filter Housing with Test Port - 10" (Req 1 - RO Carbon Block Cartridge 10" (MROTCB10-2.5-10))	
MFH20	Carbon Pre-Filter Housing - 20" (Req 1 - RO Carbon Block Cartridge 20" (MROTCB20-2.5-10))	
MFH20D	Carbon Pre-Filter Double Housing - 20" (Req 2 - RO Carbon Block Cartridges 20" (MROTCB20-2.5-10))	
MFH20DTP	Carbon Pre-Filter Double Housing with Test Port - 20" (Req 2 - RO Carbon Block Cartridges 20" (MROTCB20-2.5-10))	
MFH20T	Carbon Pre-Filter Triple Housing - 20" Req 2 - RO Carbon Block Cartridges 20" (MROTCB20-2.5-10) and 1 - RO Sediment Cartridge 20" (MROTSC20-2.5-5)	
MROPICC	Filter Element - 10" Inline Carbon Cartridge (for Education Rack)	
MROTCB10-2.5-10	Filter Element - 10" Carbon Block, 2.5" Dia, 10 Micron (for HP 150, MRO-0040, MRO-0080, MROP-0050, MFH10, MFH10D, MFH10T)	
MROTSC10-2.5-5	Filter Element - 10" Sediment, 2.5" Dia, 5 Micron (for HP 150, MRO-0040, MRO-0080, MROP-0050, MFH10, MFH10D, MFH10T)	

Water Makers



MROT-0500

Part No.#	Description	Price
MROP-0050	RO Water Maker, PuroTech 50 GPD with Dual Prefilter	
MROT-0200	RO Water Maker, Flexeon 200 GPD 115/220 Volt, 50/60 Hertz with 10" x 2.5" Triple Prefilter	
MROT-0500	RO Water Maker, Flexeon 500 GPD 115/220 Volt, 50/60 Hertz with 10" x 4.5" Dual Prefilter	
MROT-1000	RO Water Maker, Flexeon 1000 GPD 115/220 Volt, 50/60 Hertz with 10" x 4.5" Dual Prefilter	
MROT-1500	RO Water Maker, Flexeon 1500 GPD 115/220 Volt, 50/60 Hertz with 20" x 4.5" Dual Prefilter	
MROT-1800	RO Water Maker, Flexeon 1800 GPD 115/220 Volt, 50/60 Hertz with 20" x 4.5" Dual Prefilter	
MROT-2000	RO Water Maker, Flexeon 2000 GPD 115/220 Volt, 50/60 Hertz with 20" x 4.5" Dual Prefilter	
MROT-4000	RO Water Maker, Flexeon 4000 GPD 220V, 60 Hertz Pump with 20" x 4.5" Triple Prefilter (2 Sediment, 1 Carbon)	

Water Maker Replacement Membranes & Parts




Part No.#	Description	Price
MROEPK	RO End Plug Kit	
MROFCW	Canister Wrench for 10" and 20" Filter Canister	
MROCSW	Canister Strap Wrench for All Filters	
MROFSV	Feed Solenoid Valve	
MROLPS	Low Pressure Switch	
MROMTFC	Membrane, TFC HF for MRO-0025	
MROM-0040	Membrane for MRO-0040	
MROM-0300	Membrane for MRO-0300	
MROM-0600	Membrane for MRO-0600, MRO-0900 & MRO-1200 (MRO-0900 & MRO-1200 Requires 2)	
MROPM-0050	Membrane for MROP-0050 1.8" x 12"	
MROTM-0200	Membrane for MROT-0200 2.5" x 14"	
MROTM-0500	Membrane for MROT-0500 (Requires 2) and MROT-1000 (Requires 3) 2.5" x 21"	
MROTM-1500	Membrane for MROT-1500 2.5" x 40" (Requires 2)	
MROA205803	Membrane Snap Ring, 2.5"	
MROA205804	Membrane Snap Ring, 4.0"	
MROP600	Pump for HP 150-600 (Brass-Std) Does not include motor	
MROP900	Pump for HP 900/1200 (Brass-Std) Does not include motor	
MROTP-0500	Pump: Vane, 2.3 GPM for MROT-0500 RO Water Maker Does not include motor	
MROTP-1000	Pump: Vane, 3 GPM for MROT-1000 & MROT-1500 RO Water Maker Does not include motor	

Water Maker Replacement Membranes & Parts Continued

Part No.#	Description	Price
MROTP-1800	Pump: Vane, 5.3 GPM for MROT-1800 & MROT-2000 RO Water Maker	
MROPM033	Motor 1/3 HP for HP 150-600	
MROPM050	Motor 1/2 HP for HP 900-1200	
MROTPM-0500	Motor 1/3 HP for MROT-0500	
MGSSBA300-025	S.S. Pressure Gauge 0-300 PSI 2.5" Dia	
MGSS15-025	S.S. Pressure Gauge 0-15 PSI 2.5" Dia	
MGSS30-025	S.S. Pressure Gauge 0-30 PSI 2.5" Dia	
MROSK110	Solenoid Kit Includes: Solenoid Valve, 1" FPT, DIN and 3 Meter Cord. 110V. (For MROT-0500, MROT-1000, MROT-1500, MROT-1800, MROT-2000)	
MROSK220	Solenoid Kit Includes: Solenoid Valve, 1" FPT, DIN and 3 Meter Cord. 220V. (For MROT-4000)	
MROSK-0200	Solenoid Valve, 3/8" BSP, 24V (For MROT-0200)	

**Membranes priced individually.

RO Filters (Replaceable Elements)

Part No.#	Description	Price
MROSCP	Filter Element - 10" Sediment/Carbon Prefilter (for MRO-0025)	
MROBCF	Filter Element - 10" Carbon Postfilter (for MRO-0025)	
MROTCB10-2.5-10	Filter Element - 10" Carbon Block, 2.5" Dia, 10 Micron	
MROTCB10-4.5-10	Filter Element - 10" Carbon Block, 4.5" Dia, 10 Micron 	
MROTCB20-2.5-10	Filter Element - 20" Carbon Block, 2.5" Dia, 10 Micron	
MROTCB20-4.5-10	Filter Element - 20" Carbon Block, 4.5" Dia, 10 Micron	
MROTSC10-2.5-5	Filter Element - 10" Sediment, 2.5" Dia, 5 Micron	
MROTSC10-4.5-5	Filter Element - 10" Sediment, 4.5" Dia, 5 Micron 	
MROTSC20-2.5-5	Filter Element - 20" Sediment, 2.5" Dia, 5 Micron	
MROTSC20-4.5-1	Filter Element - 20" Sediment, 4.5" Dia, 1 Micron	
MROTSC20-4.5-5	Filter Element - 20" Sediment, 4.5" Dia, 5 Micro 	
MRO-TFHO	O-Ring for Triple Filter Housing	

Dosing

Dosing & Control

Part No.#	Description	Price
DD03	Chem Feed Dosing & Control For pH and Conductivity with Controllers, Pumps, (2) 3 Gallon Tanks. Rack or Storage Tank Mounted. (up to 4 racks)	
DD05	Chem Feed Dosing & Control For pH and Conductivity with Controllers, Pumps, (2) 5 Gallon Barrels and Stainless Steel Stand. (3-8 racks)	
DD15	Chem Feed Dosing & Control For pH and Conductivity with Controllers, Pumps, (2) 15 Gallon Barrels and Stainless Steel Stand. (12-24 racks)	
DD30	Chem Feed Dosing & Control For pH and Conductivity with Controllers, Pumps, (2) 30 Gallon Barrels and Stainless Steel Stand. (25-50 racks)	

Stacked Dosing & Control

Part No.#	Description	Price
DDS05	Chem Feed Dosing & Control For pH and Conductivity with Controllers, Pumps, (2) 5 Gallon Stacked Barrels and Stainless Steel Stand. (3-8 racks)	
DDS15	Chem Feed Dosing & Control For pH and Conductivity with Controllers, Pumps, (2) 15 Gallon Stacked Barrels and Stainless Steel Stand. (12-24 racks)	

Barrels & Pumps

Part No.#	Description	Price
DAD03	Chem Feed Dosing Tanks & Pumps For pH and Conductivity with (2) 3 Gallon Tanks (Requires Controllers)	
DAD05	Chem Feed Dosing Barrels & Pumps For pH and Conductivity with (2) 5 Gallon Barrels (Requires Controllers)	
DAD15	Chem Feed Dosing Barrels & Pumps For pH and Conductivity with (2) 15 Gallon Barrels (Requires Controllers)	
DAD30	Chem Feed Dosing Barrels & Pumps For pH and Conductivity with (2) 30 Gallon Barrels (Requires Controllers)	

Barrels & Pumps

Part No.#	Description	Price
DAS05	Chem Feed Dosing Barrel & Pump For pH and Conductivity with (1) 5 Gallon Barrel (Requires Controllers)	
DAS15	Chem Feed Dosing Barrel & Pump For pH and Conductivity with (1) 15 Gallon Barrel (Requires Controllers)	
DAS30	Chem Feed Dosing Barrel & Pump For pH and Conductivity with (1) 30 Gallon Barrel (Requires Controllers)	

Stacked Barrels & Pumps

Part No.#	Description	Price
DADS05	Chem Feed Dosing Stacked Barrels & Pumps For pH and Conductivity with (2) 5 Gallon Barrels (Requires Controllers)	
DADS15	Chem Feed Dosing Stacked Barrels & Pumps For pH and Conductivity with (2) 15 Gallon Barrels (Requires Controllers)	
DADS30	Chem Feed Dosing Stacked Barrels & Pumps For pH and Conductivity with (2) 30 Gallon Barrels (Requires Controllers)	

Controllers

Part No.#	Description	Price
MCAPhCond	Complete Control Assembly: Signet Universal Controller (for pH and Conductivity) including Probes, NEMA 4 Enclosure, Wiring (To dose pH and Salinity, requires pumps & barrels)	
MCAPhCond-S	Complete Control Assembly for Seawater: Signet Universal Controller (for pH and Conductivity - Seawater Ranges 100us to 200,000us) including Probe, NEMA 4 Enclosure, Wiring (To dose pH and Salinity, requires pumps & barrels)	

Part No.#	Description	Price
MCHpHCdHH	Handheld pH / Conductivity Meter	
MCHpHCD-Probe	Replacement Electrode for Hanna Hand Held pH/Conductivity Meter	
MCGF9900	Universal Controller for pH or Conductivity (3-9900-1P)	
MCG3-9900.394	Direct Conductivity/Resistivity Module for Universal Signet Controller (3-900.94)	
MCG3-9900.393	Relay Module for Universal Signet Controller (3-900.93)	
MCG3-8050.396	RC Filter Kit for Universal Signet Controller (3-8500.396)	
MCPAPH-2750-3	pH Pre-Amp (Submersible Dry Lock) FOR SYSTEMS WITH SIGNET 3-9900-1P CONTROLLER:	
MCPAPH-2760	pH Pre-Amp (Dry Lock) FOR SYSTEMS WITH SIGNET 8750-2P CONTROLLER:	
MCEPH-2724	pH Electrode, Black (Dry Lock)	
MCECOND	Conductivity Electrode 1.0 cell for 0-15 ppt	
MCSPH7	Calibration Solution pH 7 Packet	
MCSPH4	Calibration Solution pH 4 Packet	
MCSPH7C	Calibration Solution pH 7 Case	
MCSPH4C	Calibration Solution pH 4 Case	
MCSC0718	Calibration Solution Conductivity Bottle for Aquadyne	
MCSpHAQD	Calibration Solution pH and ORP for Aquadyne	
MCSC1413	Calibration Solution Conductivity Packet	
MCSC1413C	Calibration Solution Conductivity Case	

Storage Tanks

Part No.#	Description	Price
MSTRCT015	15 Gallon Storage Tank w/Stand, Float Valve and Strut to Mount DD03 and/or Small RO	
MSTRND030	30 Gallon Plastic Storage Tank w/ Stand & Float Valve, Connection to Sump Tank. 20" Diameter x 30" Tall (46" Tall with Stand)	
MSTRCT030	30 Gallon Storage Tank with Stand & Float Valve, Connection to Sump Tank. Square Plastic 12" x 12" x 48" Tall (64" Tall with Stand)	
MSTRND050	50 Gallon Plastic Storage Tank w/ Stand & Float Valve, Connection to Sump Tank. 22" Diameter x 34" Tall (50" Tall with Stand)	
MSTRCT075	75 Gallon (284 Liter) Plastic Storage Tank w/ Stand & Float Valve, Connection to Sump Tank 18" x 18" x 60" Tall (76" Tall with Stand)	
MSTRND100	100 Gallon Plastic Storage Tank w/ Stand & Float Valve, Connection to Sump Tank. Round Plastic. 27" Diameter x 48" Tall (64" Tall with Stand)	
MSTRCT140	140 Gallon Storage Tank w/ Stand & Float Valve, Connection to Sump Tank. Square Plastic. 24" x 24" x 58" Tall (74" with Stand)	
MSTRCT200	200 Gallon Storage Tank w/ Stand & Float Valve, Rectangular Aluminum Powder Coated Construction. 24" x 48" x 48" Tall	
MSTRCT300	300 Gallon Storage Tank w/ Stand and Float Valve, Rectangular Aluminum Powder Coated Construction. 48" x 24" x 70" Tall.	
MSTRND400	400 Gallon Plastic Storage Tank w/ Stand and Float Valve, Connection to Sump Tank. Round Plastic. 39 Diameter" x 84" Tall (96" Tall with Stand)	
MSTRCT005-SP	Stand Alone Rack 5 Gallon Overflow Tank with Automatic Sump Pump (Danner 1200) to Connect to Sink.	

Note: Custom Sizes Available

Timed Water Exchange

Part No.#	Description	Price
MCTWX05	Automatic Timed Water Exchange 5 GPM with Danner 350 Pump with MCTWEDRI Timer Box, Anti-siphon Valve and Hose	
MCTWX05EU	Automatic Timed Water Exchange 20 LPM with Quiet One 3000 Pump with MCTWEDRI Timer Box, Anti-siphon Valve and Hose. (SPECIFY COUNTRY & PLUG TYPE)	
MCTWX10	Automatic Timed Water Exchange 10 GPM with Danner 1200 Pump with MCTWEDRI Timer Box, Anti-siphon Valve and Hose	
MCTWX10EU	Automatic Timed Water Exchange 40 LPM with Quiet One 4000 Pump with MCTWEDRI Timer Box, Anti-siphon Valve and Hose. (SPECIFY COUNTRY & PLUG TYPE)	
MCTWX25	Automatic Timed Water Exchange 25 GPM with March TE5.5 Pump with MCTWEDRI Timer Mounted in Box on Tank (Pump Stand for Large Aluminum Storage Tanks Only)	
MCTWX25EU	Automatic Timed Water Exchange 100 LPM with March TE5.5 Pump with MCTWEDRI Timer Mounted in Box on Tank (Pump Stand for large aluminum storage tanks only.) (SPECIFY COUNTRY & PLUG TYPE)	
MCTWXD1800	Automatic Timed Water Exchanged 15 GPM with Danner 1800 Pump with MCTWEDRI Timer Box, Anti-siphon Valve and Hose	
MCTWXD2600	Automatic Timed Water Exchanged 25 GPM with Danner 2600 Pump with MCTWEDRI Timer Box, Anti-siphon Valve and Hose	
MSTRPCDT12	Packed Column Degassing Tower 12 GPM with Danner 1800 Pump and 12" x 12" x 24" Packed Column and Mounting Hardware	
MCTWEDRCD-120VAC	Push-Button Digital Count Down Timer for Automatic Water Exchange 120 Volt Pump	
MCTWEDRG-120VAC	Timer Control Box with Grasslin Timer, Outlet and Cover 120VAC	
MCTWEDRG-240VAC	Timer Control Box with Grasslin Timer, Outlet and Cover 240VAC	
MCTWEDRI-120VAC	Timer Control Box with Intermatic Timer, Outlet and Cover 120VAC	
MCTWEVSG-120VAC	Timer Control Box for Solenoid Valve with Grasslin Timer 120VAC with 24VAC Output	
MCTWEVSG-240VAC	Timer Control Box for Solenoid Valve with Grasslin Timer 240VAC with 24VAC Output	
MCTLG-120VAC	Timer Control Box for Light Box with Grasslin Timer 120VAC	
MCTLG-240VAC	Timer Control Box for Light Box with Grasslin Timer 240VAC	
ETG-120VAC	Intermatic/Grasslin Timer - 120VAC Replacement Timer Switch for MCTWEDRG-120VAC	
ETG-240VAC	Intermatic/Grasslin Timer - 240VAC Replacement Timer Switch for MCTWEDRG-240VAC	
ETIDT620	Plug-in Timer - 120VAC	
IPDC-2S	Rack to Rack Connection Package for 2 Single Sided Racks. Plumbing and connections includes 1" Flexible Hose, bulkhead fittings, 3/8" Polyethylene Tube and Danner 250 GPH Pump. Note: Bulkheads pre-installed at factory	
IPDC-2SEU	Rack to Rack Connection Package for 2 Single Sided Racks. Plumbing and connections includes 25mm Flexible Hose, bulkhead fittings, 10mm Polyethylene Tube and Quiet One 3000 Pump. Note: Bulkheads pre-installed at factory (SPECIFY COUNTRY & PLUG TYPE)	
ECPS1002715	Pump Float Switch Pump Down 15' Cord w/Plug 120V	
ECPS1002721	Pump Float Switch Pump Up 15' Cord w/Plug 120V	

Monitoring and Alarms with Dosing Control

Part No.#	Description	Price
MMCAQT5000	<p>Aquatouch Monitoring System with pH and Conductivity Controllers and NEMA 4 Enclosure</p> <p>Monitors: pH, Conductivity, Temperature, Water Level and Filter Pressures. (To dose pH and salinity, requires barrels and pumps)</p> <p>Aquatouch features:</p> <p>Controls: Turns off System Pumps and UV in Case of Low Water Level</p> <p>Alarms: SMS and E-mail Alarms Port Forwarding Needs to be Established with Each Monitoring Unit Given a MAC Address to Allow the Aquatics Monitoring System to Send Out Email and Text Message Alerts</p> <p>Security: Multi-Level Security Built in Allows Users to See Parameters and Only Administrators to Make Changes</p> <p>Remote Access: Built in Web Server Allows Programming Updates and Changes from any Computer with Internet Access. Multi-Level Alarm Capability</p> <p>Data-log: Micro SD Card (up to 2 GB) Logs Data for Download via Web Server</p>	
MMCS3000	<p>Sensaphone SCADA 3000 Monitoring System with pH and Conductivity Controllers and NEMA 4 Enclosure</p> <p>Monitors: pH, Conductivity, Temperature, Water Level and Filter Pressures. (To dose pH and salinity, requires barrels and pumps)</p> <p>SCADA 3000 features:</p> <p>Multi-Level Phone, Pager, Fax and E-mail Warnings and Alarms</p> <p>Pre-Programmed User-Defined Settings</p> <p>Real-Time Remote Monitoring Software</p> <p>Digital Display with LED Indicator Lights</p> <p>Can Be Connected to BAS</p> <p>Web Server Remote Login for Viewing System Parameters in Real Time</p> <p>Up to 64 Alarm Recipients and 16 Input Capacity (requires dedicated analog phone line)</p>	
MMCW600	<p>Web 600 Monitoring System with pH and Conductivity Controllers</p> <p>Monitors: pH, Conductivity, Temperature and Water Level. (To dose pH and salinity, requires barrels and pumps)</p> <p>Sensaphone Model Web 600 features:</p> <p>Data-Logging</p> <p>SMS (Text Messages) and E-mail Alarms</p> <p>Web Server Remote Login for Viewing System Parameters in Real Time</p> <p>Up to 8 Alarm Recipients and 6 Input Capacity (requires Static IP Address and RG45 Ethernet Conn)</p>	
MMAQT5000	<p>Monitoring System, Aquatouch</p> <p>Monitors: pH, Conductivity, Water Temperature, Room Temperature and Low Water Level.</p> <p>Add-ons: Filter Pressure and Dissolved Oxygen.</p> <p>Controls: Turns off System Pumps and UV in case of Low Water Level.</p> <p>Alarms: SMS and E-mail alarms. Port Forwarding Needs to be Established with Each Monitoring Unit Given a MAC Address to Allow the Aquatics Monitoring System to Send Out Email and Text Message Alerts</p> <p>Security: Multi-Level Security built in allows users to see parameters and only administrators to make changes.</p> <p>Remote Access: Built-in Web Server allows programming updates and changes from any computer with Internet access. Multi-Level Alarm capability.</p> <p>Data-log: Micro SD card (up to 2 GB) logs data for download via Web Server.</p>	



Part No.#	Description	Price
MMS3000	<p>Monitoring and Alarm System w/ NEMA 4 Enclosure Sensaphone SCADA 3000 Pre-programmed User-defined Settings Real-Time Remote Monitoring Software Digital Display with LED Indicator Lights Can be connected to BAS Monitors pH, Temperature, Conductivity, and Water Level Pressure Sensors on Finish Filters Multi-Level Phone, Pager, Fax and Email Warnings and Alarms Web Server Remote Login for viewing system parameters in real time Up to 64 Alarm Recipients and 16 Input Capacity (requires dedicated analog phone line)</p>	
MMS0600	<p>Monitoring System Monitors pH, Temperature, Conductivity, and Water Level (pH and conductivity only when connected to Aquaneering Dosing System) Sensaphone Model Web 600 features: Data-Logging SMS (Text Messages) and E-mail Alarms Web Server Remote Login for viewing system parameters in real time Up to 8 Alarm Recipients, and 6 Input Capacity Requires Static IP Address and RG45 Ethernet connection Port Forwarding Needs to be Established with Each Monitoring Unit Given a MAC Address to Allow the Aquatics Monitoring System to Send Out Email and Text Message Alerts</p>	

Monitoring Accessories

Part No.#	Description	Price
MMS3000RT	Room Temperature Sensor	
MMS3000WT	Water Temperature Sensor	
MMS3000BAT	Back-Up Battery for SCADA 3000	
MMS3000DO	Dissolved Oxygen Sensor	
MMS-3000DOMK	Dissolved Oxygen Sensor Membrane Kit for 5 membrane changes and solution	
MMS3000DO-SPK	Dissolved Oxygen Spare Parts Kit includes 2 membrane tensioning washers, 2 membrane locks and 1 membrane replacement tool	
MMS3000PT	Pressure Transmitter Sensor (0-50 psi)	
MMS3000SL	Sump Level Sensor	
MMS3000EM	Expansion Module	
MMS3000TGP	Total Gas Pressure Probe	
MMS3000T-	Total Gas Pressure Probe Membrane Cartridge	
MMS3000TG-	Total Gas Pressure Cartridge Replacement O-Ring	
MMS3000PS	Power Supply 115/240 , 50-60 Hz for SCADA 3000	
MMS3000EK	Ethernet Kit for SCADA 3000	
MMS3000WD	Water Detection Sensor for Floor to connect to MAQT5000, MS3000, MS1800 and MS0600	

Ultraviolet Sterilizer

Part No.#	Description	Price
MUV008	8 Watt UV Sterilizer 120 V 50/60 Hz Includes 8 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV018	18 Watt Smart UV Sterilizer 120 V 50/60 Hz Includes 18 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV025	25 Watt Smart UV Sterilizer 120-240 V ac 50/60 Hz Includes 25 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV025E	25 Watt Smart UV Sterilizer for Europe 230 V 50/60 Hz Includes 25 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV040	40 Watt Smart UV Sterilizer 120-240 V ac 50/60 Hz Includes 40 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV040E	40 Watt Smart UV Sterilizer for Europe 230 V 50/60 Hz Includes 40 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV080	80 Watt Smart High Output UV Sterilizer 120 V 50/60 Hz Includes 80 Watt UV Lamp, Quartz Sleeve and Ballast	
MUV080E	80 Watt Smart High Output UV Sterilizer for Europe 230 V 50/60 Hz Includes 80 Watt UV Lamp, Quartz Sleeve and Ballast	

Note: Add "-B" to end of part number to designate bayonet plug. Add additional \$28 to price.

Commercial UV Sterilizer

Aquaneering Commercial 100-200 Watt UV Sterilizers use the latest in Amalgam Lamp Technology.

Part No.#	Description	Price
MUVAQ100-*	100 Watt UV Aquaneering Commercial UV 3" Body x 40" L; 2" Inlet/Outlet Flow Rate: 25 GPM 120V, 230V 50/60Hz Includes 100 Watt UV Lamp, Quartz Sleeve and Ballast	
MUVAQ200-*	200 Watt UV Aquaneering Commercial UV 4" Body x 52" L; 2" Inlet/Outlet Flow Rate: 50 GPM 120V, 230V 50/60Hz Includes 200 Watt UV Lamp, Quartz Sleeve and Ballast	
MUVAQ400-*	400 Watt UV Aquaneering Commercial UV 4" Body Dbl x 52" L; 3" Inlet/Outlet Flow Rate: 100 GPM 120V, 230V 50/60Hz Includes Two 200 Watt UV Lamps, Quartz Sleeves and Ballasts	
MUVAQ600-*	600 Watt UV Aquaneering Commercial UV 4" Body Triple x 52" L; 4" Inlet/Outlet Flow Rate: 200 GPM 120V, 230V 50/60Hz Includes Three 200 Watt UV Lamps, Quartz Sleeves and Ballasts	

Note: Add "-120" or "-230" to end of part number to designate voltage. Add additional \$30 for 230v.

Stainless Steel UV Sterilizer

Part No.#	Description	Price
MUVSS025	25 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 2-3/8" Dia, 1" Inlet/Outlet, 25" Long 120V - 277V 50/60Hz Includes 25 Watt UV Lamp, Quartz Sleeve and Ballast	
MUVSS040	40 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 2-3/8" Dia, 1" Inlet/Outlet, 40" Long 120V - 277V 50/60Hz Includes 40 Watt UV Lamp, Quartz Sleeve and Ballast	
MUVSS080	80 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 3-1/2" Dia, 1" Inlet/Outlet, 40" Long 120V - 277V 50/60Hz Includes 80 Watt UV Lamp, Quartz Sleeve and Ballast	

Note: Add "-B" to end of part number to designate bayonet plug. Add additional \$28 to price.

Part No.#	Description	Price
MUVSS100-*	100 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 3-1/2" Dia, 1-1/2" Inlet/Outlet, 38" Long 120V or 230V, 50/60Hz Includes 100 Watt UV Lamp, Quartz Sleeve and Ballast	
MUVSS200-*	200 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 4-1/2" Dia, 2" Inlet/Outlet, 50" Long 120V or 230V, 50/60Hz Includes 200 Watt UV Lamp, Quartz Sleeve and Ballast	
MUVSS400-*	400 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 8-5/8" Dia, 3" Inlet/Outlet, 50" Long 120V or 230V, 50/60Hz Includes Two 200 Watt UV Lamps, Quartz Sleeves and Ballasts	
MUVSS600-*	600 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 8-5/8" Dia, 2" Inlet/Outlet, 50" Long 120V or 230V, 50/60Hz Includes Three 200 Watt UV Lamps, Quartz Sleeves and Ballasts	
MUVSS800-*	800 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 8-5/8" Dia, 4" Inlet/Outlet, 50" Long 120V or 230V, 50/60Hz Includes Four 200 Watt UV Lamps, Quartz Sleeves and Ballasts	
MUVSS1000-*	1000 Watt UV Sterilizer Stainless Steel Manufactured by Aquaneering 8-5/8" Dia, 4" Inlet/Outlet, 50" Long 120V or 230V, 50/60Hz Includes Five 200 Watt UV Lamps, Quartz Sleeves and Ballasts	

Note: Add "-120" or "-230" to end of part number to designate voltage. Add additional \$30 for 230v.

UV Replacement Parts

UV Replacement Housing

Part No.#	Description	Price
MUVAQ100H	100 Watt Replacement Housing for Aqua UV (does not include lamp, sleeve or ballast box)	
MUVAQ200H	200 Watt Replacement Housing for Aqua UV (does not include lamp, sleeve or ballast box)	
MUVAQ400H	400 Watt Replacement Housing for Aqua UV (does not include lamp, sleeve or ballast box)	
MUVAQ600H	600 Watt Replacement Housing for Aqua UV (does not include lamp, sleeve or ballast box)	
MUVAQ800H	800 Watt Replacement Housing for Aqua UV (does not include lamp, sleeve or ballast box)	
MUVSSOR	Aquaneering Stainless Steel UV Sterilizer O-Ring	

UV Replacement Lamps

Part No.#	Description	Price
MUL008AU	8 Watt UV Replacement Lamp for Aqua UV	
MUL015AU	15 Watt UV Replacement Lamp for Aqua UV (11.31"L)	
MUL025	25 Watt UV Replacement Lamp for Emperor, Aqua UV and Aquaneering (17.25"L)	
MUL030	30 Watt UV Replacement Lamp, Aquanetics (35.13"L)	
MUL040	40 Watt UV Replacement Lamp for Emperor, Aqua UV, and Aquaneering (33.25" L)	
MUL057	57 Watt UV Replacement Lamp for Emperor (18.25"L)	
MUL057AU	57 Watt UV Replacement Lamp for Aqua UV (17.19"L)	
MUL065	65 Watt UV Replacement Lamp for Emperor (61.75"L)	
MUL080	80 Watt UV Replacement Lamp for Emperor and Aquaneering (38.25" L)	
MUL100	100 Watt UV Replacement Lamp for Aqua UV and Aquaneering Commercial (36"L)	
MUL150	150 Watt T6 Style UV Replacement Lamp for Emporer (60.75"L)	
MUL200	200 Watt UV Replacement Lamp for Aqua UV and Aquaneering Commercial (48"L)	

UV Replacement Ballasts

Part No.#	Description	Price
MUBA040AQ	UV Ballast Assembly: 25, 40 Watt in Enclosure 120-240v, 50/60 Hz	
MUBA040AQ-B	UV Ballast Assembly: 25, 40 Watt in Enclosure 120-240v, 50/60 Hz with Bayonet Plug	
MUBA080AQ	UV Ballast Assembly: 80 Watt in Enclosure 120-240v, 50/60 Hz	
MUBA080AQ-B	UV Ballast Assembly: 80 Watt in Enclosure 120-240v, 50/60 Hz with Bayonet Plug	
MUBA100AQ-120 *	UV Ballast Assembly: 100-200 Watt in Enclosure 120v, 60 Hz	
MUBA100AQ-230	UV Ballast Assembly: 100-200 Watt in Enclosure 230v, 50/60 Hz	
MUBA400AQ-120 *	UV Ballast Assembly: 400 Watt (For 2 200 Watt Lamps) in Metal Gutter 120v	
MUBA400AQ-230	UV Ballast Assembly: 400 Watt (For 2 200 Watt Lamps) in Metal Gutter 230v	
MUBA600AQ-120 *	UV Ballast Assembly: 600 Watt (For 3 200 Watt Lamps) in Custom Enclosure 120v	

UV Replacement Ballasts Continued

Part No.#	Description	Price
MUBA600AQ-230	UV Ballast Assembly: 600 (For 3 200 Watt Lamps) in Custom Enclosure 230v	
MUBA800AQ-120 *	UV Ballast Assembly: 800 Watt (For 4 200 Watt Lamps) in Custom Enclosure 120v	
MUBA800AQ-230	UV Ballast Assembly: 800 Watt (For 4 200 Watt Lamps) in Custom Enclosure 230v	
MUBA1000AQ-120 *	UV Ballast Assembly: 1,000 Watt (For 5 200 Watt Lamps) in Custom Enclosure 120v	
MUBA1000AQ-230	UV Ballast Assembly: 1,000 Watt (For 5 200 Watt Lamps) in Custom Enclosure 230v	
MUBA1200AQ-120 *	UV Ballast Assembly: 1,200 Watt (For 6 200 Watt Lamps) in Custom Enclosure 120v	
MUBA1200AQ-230	UV Ballast Assembly: 1,200 Watt (For 6 200 Watt Lamps) in Custom Enclosure 230v	

Note: For 40 and 80 Watt ballasts, add "-B" to end of part number to designate bayonet plug. Add additional \$28 to price.

UV Replacement Quartz Sleeve

Part No.#	Description	Price
MUQ008AU	8 Watt UV Replacement Quartz Sleeve White Aqua UV	
MUQ025	25 Watt UV Replacement Quartz Sleeve for Emperor (23.5"L)	
MUQ030	30 Watt UV Replacement Quartz Sleeve for Aquanetics (35.13"L)	
MUQ040	40 Watt and 80 Watt UV Replacement Quartz Sleeve for Emperor (38.25"L)	
MUQ065	65 Watt and 150 Watt UV Replacement Quartz Sleeve for Emperor (65"L)	
MUQ015AU	15 Watt UV Replacement Quartz Sleeve White Aqua UV and Aquaneering Commercial UV (15"L)	
MUQ025AU	25 Watt UV Replacement Quartz Sleeve White Aqua UV (19"L)	
MUQ040AU	40 Watt UV Replacement Quartz Sleeve White Aqua UV (38.25"L)	
MUQ100	100 Watt UV Replacement Quartz Sleeve (36"L)	
MUQ200	200 Watt UV Replacement Quartz Sleeve (48"L)	
MUSK20375	Service Kit for 18, 25 and 40 Watt UV Sterilizers Includes Power Supply Gasket, O-Ring for Quartz Sleeve, Retaining Module, and Quartz Sleeve Gasket	
MUSK20375-2	Service Kit for 80 Watt Smart HO UV Sterilizers Includes Power Supply Gasket, O-Ring for Quartz Sleeve, Retaining Module, and Quartz Sleeve Gasket	
MUSK-COMCL	Service Kit for 100/150/200 Watt Smart HO UV Commercial Sterilizer (including post 2009 multiple lamp ((ML)) units) Includes Quartz Sleeve O-Ring and Quartz Sleeve Gasket	
MUSK030	30 Watt O-Ring Service Kit (Aquanetics)	

Fiber Optic Parts

Part No.#	Description	Price
MPG2000L150	Replacement Lamp for Fiber Optic Light Source	

Acrylic Lighting

Part No.#	Description	Price
ECPLWALA32	32 Watt, 52" Acrylic Light (includes tube and lamp)	
ECPWALA40	40 Watt, 66" Acrylic Light (includes tube and lamp)	
ECPLWALA24	24 Watt, 17" Acrylic Light (includes tube and lamp)	
ECPLBAWM32	32 Watt Replacement Ballast for Acrylic Light	
ECPLBAWM24	24 Watt Replacement Ballast for Acrylic Light	
ECPLWM-LA32T8	32 Watt Lamp for Acrylic Light	
ECPLWM-LA24T8	24 Watt Lamp for Acrylic Light	
ECPLWALA32-T	32 Watt, 52" Acrylic Tube	
ECPLWALA24-T	24 Watt, 17" Acrylic Tube	

LED Lighting

Part No.#	Description	Price
ECPLEDRP120	LED Light Rope 120V 60Hz	
ECPLEDRPPC	LED Light Rope Power Cord	
ECPLEDRPCC	LED Light Rope Cord Connector	
ECPLEDCLED-12V66RL-WW	LED Light Rope for Light Enclosures 12v, DC, "Warm White" (74' enough for 5 or 6 shelves)	
ECPLEDCLED-12V66RL-WW-TT	LED Light Rope for Light Enclosures 12v, DC, "Warm White" (37' enough for 3 shelves)	
ECZLH220	220V Power Source and Timer for LED Light Rope	
ECZLH115	115V Power Source and Timer for LED Light Rope	
ECPLEDS-BST24-024	LED Light Strip, 24" Length, 24V	
ECPLEDBL12	Replacement 3" Flood Lamp for Zebrafish Spawning Chamber	

System Chillers

Part No.#	Description	Price
MALCHY14	Chiller with Controller, Drop In Titanium Coil, 1/4 HP 115v / 60 Hz	
MALCHY13	Chiller with Controller, Drop In, Titanium Coil, 1/3 HP 115v / 60 Hz	
MALCHY12	Chiller with Controller, Drop In Titanium Coil, 1/2 HP 115v / 60 Hz	
MALCHY12EU	Chiller with Controller, Drop In Titanium Coil, 1/2 HP 230v / 50 Hz with European Plug	
MALCHY12UK	Chiller with Controller, Drop In Titanium Coil, 1/2 HP 230v / 50 Hz with UK Plug	
MALCHDS7	Chiller, Titanium Tube In Shell, Air Cooled, 3/4 HP 115v / 60 Hz	
MALCHDS8	Chiller, Titanium Tube In Shell, Air Cooled, 3/4 HP 230v / 50 Hz	
MALCHDS9	Chiller, Titanium Tube In Shell, Air Cooled, 1 HP 230v / 50 Hz	
MALCHDS9-WC	Chiller, Titanium Tube In Shell, Water Cooled, 1 HP 230v / 50 Hz	
MALCHDS10	Chiller, Titanium Tube in Shell, Air Cooled, 1.5 HP 230v / 50 Hz	
MALCHAC115	Chiller Controller for 1/6, 1/4, 1/3, 1/2 HP, 115v / 60 Hz	
MALCHAC230	Chiller Controller for 1/6, 1/4, 1/3, 1/2 HP, 230v / 50 Hz	
MCDSFS-1.0	1" Flow Switch for Delta Star Chiller up to 1/3HP (115v)	
MCDS-FS-1.5-230	1.5" Flow Switch for Delta Star Chiller 1/2HP - 1.5HP (230v)	
MCDSFS-1.5	1.5" Flow Switch for Delta Star Chiller 3/4HP (115v)	
MCSDFS-015	1.5" Flow Switch for Delta Star Chiller 1/2HP (115v)	

System Heaters

Part No.#	Description	Price
MH100GC	Heater 100 Watt Drop In, Glass, with Built-In Controller for 20L Brine Shrimp Hatcher	

Titanium Heaters

Part No.#	Description	Price
MHT300	300 Watt Titanium Heater for ZTT350	
MHT300-B	300 Watt Titanium Heater with Bayonet Plug for ZTT350	
MHT800	800 Watt Titanium Heater for ZS560	
MHT800-B	800 Watt Titanium Heater with Bayonet Plug for ZS560	
MHT800-230	800 Watt Titanium Heater for Stand Alone Rack 230v 50Hz	
MHT800-230-B	800 Watt Titanium Heater with Bayonet Plug for Stand Alone Rack, 230v 50Hz	
MHT800C	800 Watt Heater Controller 115v 60Hz	

Single Tube Bottom Heaters

Part No.#	Description	Price
MH1800L-120	1800 Watt Heater, L-Shaped, 120v (16" x 17" Tall)	
MH1800L-230	1800 Watt Heater, L-Shaped, 230v (16" x 17" Tall)	
MH2500L-208	2500 Watt Heater, L-Shaped, 208v, 1-Phase (21" x 24" Tall)	
MH2500L-230	2500 Watt Heater, L-Shaped, 230v, 1-Phase (21" x 24" Tall)	
MH3500L-208	3500 Watt Heater, L-Shaped, 208v, 1-Phase (27" x 24" Tall)	
MH3500L-230	3500 Watt Heater, L-Shaped, 230v, 50 Hz / 208v, 60 Hz 1-Phase (27" x 24" Tall)	
MH5000L-208	5000 Watt Heater, L-Shaped, 208v, 1-Phase (36" x 48" Tall)	
MH5000L-230	5000 Watt Heater, L-Shaped, 230v, 50 Hz / 208v, 60 Hz 1-Phase (36" x 48" Tall)	
MH6000L-208	6000 Watt Heater, L-Shaped, 208v, 1-Phase (44" x 48" Tall)	
MH6000L-230	6000 Watt Heater, L-Shaped, 230v, 50 Hz / 208v, 60 Hz 1-Phase (44" x 48" Tall)	

Triple Tube Bottom Heaters

Part No.#	Description	Price
MH3000L3-1P208	3000 Watt Heater, Triple Tube, L-Shaped, 240v Single-Phase (11" x 24" Tall)	
MH3000L3-3P208	3000 Watt Heater, Triple Tube, L-Shaped, 208v, 3-Phase (11" x 24" Tall)	
MH5000L3-3P208	5500 Watt Heater, Triple Tube, L-Shaped, 208v, 3-Phase (16" x 24" Tall)	
MH7500L3-3P208	7500 Watt Heater, Triple Tube, L-Shaped, 208v, 3-Phase (21" x 24" Tall)	
MH10500L3-3P208	10,500 Watt Heater, Triple Tube, L-Shaped, 208v, 3-Phase (27" x 48" Tall)	

1 1/4" NPT SCrewplug Heater

Part No.#	Description	Price
MH1000S-120	1000 Watt Heater, Screwplug, 120v, 60 Hz (10" long)	
MH1000S-120-B	1000 Watt Heater, Screwplug, with Bayonet Plug 120v, 60 Hz (10" long) for ZS and ZD Racks	
MH1000S-EZ120	1000 Watt Heater, Screwplug, EZ Plug, 120v, 60 Hz (10" long)	
MH1000S-230	1000 Watt Heater, Screwplug, 230v, 50 Hz (10" long)	
MH1000S-230-B	1000 Watt Heater, Screwplug, with Bayonet Plug 230v, 50 Hz (10" long) for ZS and ZD Racks	

Easy Plug Heater & Digital Control Systems

Part No.#	Description	Price
MHC1000-EZ120	1000 Watt Heater w/Controller, Over the Side Style, EZ Plug, 120v, 60 Hz (11"long)	
MHC1000-EZ230	1000 Watt Heater w/Controller, Over the Side Style, EZ Plug, 230v, 50 Hz (11"long)	

Easy Plug Heater [Straight / Over the Side] [No Controller]

Part No.#	Description	Price
MH1000-EZ120	1000 Watt Heater, Over the Side Style, EZ Plug, 120v, 60Hz (11" long)	
MH1000-EZ230	1000 Watt Heater, Over the Side Style, EZ Plug, 230v, 50 Hz (11"long)	

Smart One Easy Plug Metal Bottom Heater [Axial Style] [No Controller]

Part No.#	Description	Price
MH1000AXS1-EZ120	1000 Watt Heater, Axial Style, EZ Plug, 120v, 60 Hz (14" long)	
MH1000AXS1-EZ230	1000 Watt Heater, Axial Style, EZ Plug, 230v, 50 Hz (14" long)	
MH1700AXS1-EZ120	1700 Watt Heater, Axial Style, EZ Plug, 120v, 60 Hz (18" long)	
MH1700AXS1-EZ230	1700 Watt Heater, Axial Style, EZ Plug, 230v, 50 Hz (18" long)	

DRA Series Digital Thermostats [Controller]

Part No.#	Description	Price
MHCPT-DRA151	1000/1800 Watt Heater Controller, 120v, 15 Amps	
MHCPT-DRA152	1000/1800 Watt Heater Controller, 230v, 8 Amps	
MHCPT-DRA302	2500/3000/5500 Watt Heater Controller, 230 V 30 Amps, 1 or 3 phase	
MHCPT-DRA502	7500 Watt Heater Controller, 230V 50 Amps, 1 or 3 phase	
MHCPT-DRA752	10,500 Watt Heater Controller, 230 V 75 Amps, 1 or 3 phase	

Easy Plug Digital Thermostats [Controller]

Part No.#	Description	Price
MHCPTEZ120	1000/1800 Watt Heater Controller, EZ Plug, 120v, 15 Amps	
MHCPTEZ230	1000/1800 Watt Heater Controller, EZ Plug, 230v, 8 Amps	

AQ Heater Controller

Part No.#	Description	Price
MHCAQ301	Heater Controller, 120v, 30 Amps	
MHCAQ302	Heater Controller, 208v or 230v, 30 Amps	
MHCAQ502	Heater Controller, 208v or 230v, 50 Amps (For 6,000 Watt Single Phase and 10,500 Watt Three Phase)	

Parts

Part No.#	Description	Price
MFHL6021-85R	Fusible Link for L-Shaped Heaters (PL-1)	
MFHL6032-26R	Fusible Link for Over-the-Side Heaters (PT-1)	

Pumps

Pentair Pumps

Part No.#	Description	Price
MPCH050	Pump: Challenger 50 GPM ½ HP	
MPCH075	Pump: Challenger 75 GPM ¾ HP	
MPCH100	Pump: Challenger 100 GPM 1 HP	
MPCH120	Pump: Challenger 120 GPM 1.5 HP	
MPCH140	Pump: Challenger 140 GPM 2 HP	
MPCH150	Pump: Challenger 150 GPM 3 HP	
MPIF100VS	IntelliFlo VS Pump 1 HP Adjustable Variable Frequency Speed Four Programmable Flow Rates (10-90 GPM)	
MPIF300VS	IntelliFlo VS Pump 3 HP Adjustable Variable Frequency Speed Four Programmable Flow Rates (15-140 GPM)	
MPIF300VF	IntelliFlo VF Pump 3 HP Adjustable Variable Frequency Speed (15-140 GPM)	
MPWF050	Pump: WhisperFlo 1/2 HP	
MPWF075	Pump: WhisperFlo 3/4 HP	
MPWF100	Pump: WhisperFlo 1 HP	
MPWF150	Pump: WhisperFlo 1.5 HP, 208-230 Volt	
MPSH030	Pump: Seahorse 30 GPM System (MFG 38 GPM) 1/6 HP, 115V – 60 Hz	
MPSH050	Pump: Seahorse 50 GPM System (MFG 60 GPM) ½ HP, 115V – 60 Hz	
MPSH050EU	Pump: Seahorse 50 GPM System (MFG 60 GPM) ½ HP, 230V – 50 Hz	
MPSH075	Pump: Seahorse 75 GPM System (MFG 81 GPM) ¾ HP, 115V – 60 Hz	
MPSH075EU	Pump: Seahorse 75 GPM System (MFG 81 GPM) ¾ HP, 230V – 50 Hz	
MPSH100	Pump: Seahorse 100 GPM System (MFG 96 GPM) 1 HP, 115V – 60 Hz	
MPSH100EU	Pump: Seahorse 100 GPM System (MFG 96 GPM) 1 HP, 230V – 50 Hz	

Lifeguard Pumps

Part No.#	Description	Price
MPL-3000EU	Pump: QuietOne 3000 for ZT350EU and XT250EU (European Plug)	
MPL-3000EU-B	Pump: QuietOne 3000 with Bayonet Plug for ZT350EU and XT250EU	
MPL-3000PEU	Pump: QuietOne 3000 for ZT350EU and XT250EU (European Plug) (18' Chord)	
MPL-3000UK	Pump: QuietOne 3000 for ZT350UK (UK Plug)	
MPL-3000UK-B	Pump: QuietOne 3000 with Bayonet Plug for ZT350UK	
MPL-4000EU	Pump: QuietOne 4000 for ZS560EU and XS360EU Sump (European Plug) (FBB Side)	
MPL-4000EU-B	Pump: QuietOne 4000 with Bayonet Plug for ZS560EU and XS360EU Sump (FBB Side)	
MPL-4000UK	Pump: QuietOne 4000 for ZS560UK and XS360UK Sump (UK Plug)	
MPL-4000UK-B	Pump: QuietOne 4000 with Bayonet Plug for ZS560UK and XS360UK Sump	
MPL-6000EU	Pump: QuietOne 6000 for ZS560EU Rack (European Plug) (Pump Side)	
MPL-6000EU-B	Pump: QuietOne 6000 with Bayonet Plug for ZS560EU Rack (Pump Side)	
MPL-6000PEU	Pump: QuietOne 6000 for ZS560EU Rack (European Plug) (Pump Side, 18' Chord)	
MPL-6000UK	Pump: QuietOne 6000 for ZS560UK Rack (European Plug) (Pump Side)	
MPL-6000UK-B	Pump: QuietOne 6000 with Bayonet Plug for ZS560UK Rack (Pump Side)	

Lifeguard Pump Parts

Part No.#	Description	Price
MPCHPS	Pump Seal, Challenger Pump	
MPCHORL	O-Ring: Challenger, Lid	
MPLI-2200EU	Pump Impeller: Quiet One 2200/3000/4000 EU and UK	
MPLI-3000	Pump Impeller: Quiet One 3000 and 4000	
MPLI-4000HHEU	Pump Impeller: Quiet One 4000HH EU and UK	
MPLI-6000	Pump Impeller: Quiet One 4000HH and 6000	
MPLI-6000EU	Pump Impeller: Quiet One 6000 EU and UK	
MPLCS-3000	Pump Ceramic Shaft: Quiet One 3000 and 4000 (Including stabilizer top and bottom)	
MPLCS-6000	Pump Ceramic Shaft: Quiet One 4000HH and 6000 (Including stabilizer top and bottom)	
MPL-6000UK-B	Pump: QuietOne 6000 with Bayonet Plug for ZS560UK Rack (Pump Side)	

March Pumps

Part No.#	Description	Price
MPM12SB	Pump: March Magnetic Drive 12 GPM Submersible for ZD560 115V – 60 Hz	
MPM12SBEU	Pump: March Magnetic Drive 12 GPM Submersible for ZD560EU 230V – 50 Hz	
MPMSBHead	March Pump Head Assembly for Submersible Pump (Includes Impeller)	
MPMSBImp	March Pump Impeller for Submersible Pump (RH6305-655)	
MPM12AC	Pump: March Magnetic Drive 12 GPM 1/8 HP	
MPM12AC-B	Pump: March Magnetic Drive 12 GPM 1/8 HP with Bayonet Plug	
MPM12AC-RP6305.711	March Pump Head Assembly for Air Cooled Pump (Includes Impeller)	
MPM25AC	Pump: March Magnetic Drive 25 GPM 1/5 HP	
MPM30	Pump: March Magnetic Drive 30 GPM 1/3 HP	
MPM35	Pump: March Magnetic Drive 35 GPM ½ HP	

Chemical Feed Pumps

Part No.#	Description	Price
MCCFD10	Pump: Chem Feed: Dolphin 10	
MCCFD10EU	Pump: Chem Feed: Dolphin 10 230V	
MCCFS	Pump: Chem Feed: Stenner Series C45	
MCCFS45-H2	Replacement Pump Hose (2 Pieces) for Stenner Series C45 Dosing Pump	
MCCFS45-H5	Replacement Hose 5 Per Pack for C45 Stenner ChemFeed Pump	
MCCFD10PH	Replacement Head for Dolphin 10 ChemFeed Pump	
MCCFD10H	Replacement Head Tube for Dolphin 10 Chem Feed Pump	
MCCFD50	Pump: ChemFeed: Dolphin 50	
MCCFD10S	Strainer Assembly For Pulsafeeder Dolphin Pumps Includes Strainer, Coupling Connector, Cap & Ceramic Weight	
MCCFS09-H2	Replacement Pump Hose (2 Pieces) for Stenner Series 9 Dosing Pump	

Danner Pumps

Part No.#	Description	Price
MPD080	Pump: Danner Magnetic Drive 80 GPH	
MPD140	Pump: Danner Magnetic Drive 140 GPH	
MPD250	Pump: Danner Magnetic Drive 250 GPH	
MPD250Imp	Danner 250 Impeller	
MPD350	Pump: Danner Magnetic Drive 350 GPH	
MPD350HK	Danner 350 and 500 Housing Kit with O-Ring and Pump Cover	
MPD350Imp	Danner 350 and 500 Impeller	
MPD350P	Pump: Danner Magnetic Drive 350 GPH (18 ft cord)	
MPD500	Pump: Danner Magnetic Drive 500 GPH	
MPD1200	Pump: Danner Magnetic Drive 1200 GPH for ZT350 and XT250	
MPD1200-B	Pump: Danner Magnetic Drive 1200 GPH with Bayonet Plug for ZT350 and XT250	
MPD1200HK	Danner 1200 and 1800 Housing Kit with O-Ring and Pump Cover	
MPD1200Imp	Danner 1200 Impeller	
MPD1200P	Pump: Danner Magnetic Drive 1200 GPH (18 ft cord)	
MPD1800	Pump: Danner Magnetic Drive 1800 GPH for ZS560 and XT360	
MPD1800-B	Pump: Danner Magnetic Drive 1800 GPH with Bayonet Plug for ZS560 and XT360	
MPD1800Imp	Danner 1800 Impeller	
MPD1800P	Pump: Danner Magnetic Drive 1800 GPH (18 ft cord)	
MPD2400	Pump: Danner Magnetic Drive 2400 GPH for ZS660	
MPD2400-B	Pump: Danner Magnetic Drive 2400 GPH with Bayonet Plug for ZS660	
MPD2400HK	Danner 2400 and 3600 Housing Kit with O-Ring and Pump Cover	
MPD2400Imp	Danner 2400 Impeller	
MPD2600	Pump: Danner Magnetic Drive 2600 GPH for ZD Racks	
MPD2600-B	Pump: Danner Magnetic Drive 2600 GPH with Bayonet Plug for ZD Racks	
MPD3200	Pump: Danner Magnetic Drive 3200 GPH	
MPD4000	Pump: Danner Magnetic Drive 4000 GPH	
MPD4800	Pump: Danner Magnetic Drive 4800 GPH for ZD660 or larger	
MPDA20	Pump: Danner Air Pump 20 Cu In/Min	
MPDA40	Pump: Danner Air Pump 40 Cu In/Min for 20L Brine Shrimp Hatcher	
MPDA60	Pump: Danner Air Pump 60 Cu In/Min	
MPDA100	Pump: Danner Air Pump 100 Cu In/Min	
MPDSL	Pump: Danner Submersible Light with Transformer	

MDM Pumps

Part No.#	Description	Price
MPAP1221.112	Advance Pump 1/2 HP 115 Volt - 60 Hz	
MPAP1221.112EU	Advance Pump 1/2 HP 230 Volt - 50 Hz	
MPAP1341.142-XPC	Advance Pump 1 HP 115/230 Volt - 60 Hz	

Supplies

Test Kits

Part No.#	Description	Price
STK6	6-Way Test Strips (24/pack) 5-Way strip for alkalinity, pH, nitrite, nitrate, hardness. Separate ammonia strip.	
STKA	Ammonia Test Strips (24/pack)	
STKSLT	Salinity Test Strips (10/pack)	

Fish Food

Part No.#	Description	Price
SFF1	Fish Food - 1 Lb. Bag Zerbafish Select Diet (\$8.00 each for quantities 5 and up)	
SFLD0-0500	Larval Fish Food - 500g, Size 0 (20-80 micron)	
SFLD1-1000	Larval Fish Food - 1 Kilo, Size 1 (60-100 micron)	
SFBS1GRA	Artemia Brine Shrimp 90% Hatch, Grade A, 1 Lb.	
SFBS1GRB	Artemia Brine Shrimp 85% Hatch, Grade B, 1 Lb.	

Misc. Supplies

Part No.#	Description	Price
SHFNT2	2" Hagen Fish Net	
SHFNT3	3" Hagen Fish Net	
SHFNT4	4" Hagen Fish Net	
SHFNT5	5" Hagen Fish Net	
SHFNT6	6" Hagen Fish Net	
SHFN8	8" Hagen Fish Net	
SCDCST1	DeChlor-Sodium Thiosulfate 1 Pint	
SCSBC12	Sodium Bicarbonate 12 Lb.	
SCS50BA3	Sea Salt - 50 Lb. Box Crystal Sea Bioassay Formula Marine Mix	
SCAC06	Activated Carbon, Acid Washed 6 Lb. Jar 4 mm Pellets	
SCAC12	Activated Carbon, Acid Washed 12 Lb. Jar 4 mm Pellets	
ZLHSLUBE	Lubricant for Light Enclosure Door	
WPLDC111	Valve Lubricant/Sealant for FSI O-Rings 5.3 oz. Tube	

Click Pipe Clamps

Part No.#	Description	Price
PPCLK005	½" Plastic Click Pipe Clamps	
PPCLK007	¾" Plastic Click Pipe Clamps	
PPCLK010	1" Plastic Click Pipe Clamps	
PPCLK012	1 ¼" Plastic Click Pipe Clamps	
PPCLK015	1 ½" Plastic Click Pipe Clamps	
PPCLK020	2" Plastic Click Pipe Clamps	
PPCLK030	3" Plastic Click Pipe Clamps	

Bulkheads

Part No.#	Description	Price
PLB005	½" Bulkhead	
	Thread/Slip	
	Slip/Slip "SS"	
	Thread/Thread	
PLB007	¾" Bulkhead	
	Thread/Slip	
	Slip/Slip "SS"	
	Thread/ Thread	
PLB010	1" Bulkhead	
	Thread/Slip	
	Slip/Slip "SS"	
	Thread/ Thread	
PLB015	1 ½" Bulkhead	
	Thread/Slip	
	Slip/Slip "SS"	
	Thread/Thread	
PLB020	2" Bulkhead	
	Thread/Slip	
	Slip/Slip "SS"	
	Thread/Thread	

Note: Add "TS" (Thread/Slip), "SS" (Slip/Slip), or "TT" (Thread/Thread) to end of part number to designate thread/slip style.

Suction Screens & Assemblies

Part No.#	Description	Price
PLSS007-400	Plastic Pipe Larval Screen Made from PVC Pipe 3/4" - 400 Micron	
FTDLS007-400	Fish Tank Drain Tee with 2 Larval Suction Screens 3/4" - 400 Micron	
FTDASS007-LS400	Fish Tank Drain Assembly with Bulkhead, Tee, Double Suction Screens 3/4" - 400 Micron	
FTDASS007	Fish Tank Drain Assembly with Bulkhead, Overflow Tee, Suction Screen 3/4"	
FTDASS007-D	Fish Tank Drain Assembly with Bulkhead, Overflow Tee, Double Suction Screens 3/4"	

Tubing and Clamps

Part No.#	Description	Price
PTBLU025	¼" O.D. Polyurethane Blue Tube (Tank Supply)	
PTBLU025-100	¼" O.D. Polyurethane Blue Tube: 100 ft. roll	
PTBLU025-1000	¼" O.D. Polyurethane Blue Tube: 1000 ft. roll	
PTPPN025	¼" O.D. Polypropylene Natural Tube (Dosing)	
PTPPN038	3/8" O.D. Polypropylene Natural Tube	
PTAR019	3/16" O.D. Acrylic Rigid Tube: 3 ft.	
PCH038-10	3/8" Blue Self Coiling Hose 10 ft. (Nozzle Not Included)	
PCH038-25	3/8" Blue Self Coiling Hose 25 ft. (Nozzle not included)	
PVN050	1/2" Ball Valve with Female Garden Hose Connection	
PVN100	Plastic Valve Nozzle Heavy Duty Trigger 1" w/ Reducer Bushing (Hose not included)	

Contact our distributor



info@animalab.eu www.animalab.eu